EN

**ONE-DEEP WELL GUARDED PUNCH
ONE-DEEP WELL GUARDED COVER**

#17212DXXX (XXX= 750,775,800,825,850,950,1000)

#17213DYYY (YYY=950,1000)

1-INTENDED AND INDICATION FOR USE

The ONE-Deep Well Guarded Punch is specific to partially cut donor corneal tissue and is not intended for patient use. It is to be used on donor tissue only.

This system is indicated for use without restriction to the target population (age, sex and health status). There are no contra-indications

2-DESCRIPTION / FEATURES

The device is EO (Ethylene Oxide) sterilized and ready for use.

It is composed of:

- A base (A) and syringe (C): to maintain a donor cornea per aspiration (see Fig. 1)
- A blade support (B): to trephine endothelial side in different diameters (see Fig.1).

The blade support (B) with calibrated blade is available in different trephining diameters:

- Large diameters: Ø 9.50 and 10.00 mm
- Small diameters: Ø 7.50, 7.75, 8.00, 8.25 and 8.50 mm

The base (A) (where the graft is installed, centered with a groove (H) and secured to be cut) is equipped with tubing (D) to connect the syringe for aspiration.

The receptacle (see fig 2) is equipped with an aspiration groove (E) and central holes (F) to maintain the corneal graft. It also has cardinal holes (G) in case of use of a biopsy punch.

The syringe has a Luer lock connection.

The device is available in 2 versions:

- 17212Dxxx: ONE-Deep Well Guarded Punch (base, syringe and support blade in large or small diameters version)
- 17213Dxxx: ONE-Deep Well Guarded Cover (only support blade in large diameters, for additional use in combination with 17212Dxxx in small diameter)

See details of references in part "7. List of references".

3-INTENDED PERFORMANCE AND CLAIMS

More specifically, the ONE-Deep Well Guarded Punch is designed to maintain a donor cornea to:

- create one circular incision on the endothelial side at different diameters.
- manipulate a donor cornea using common hand-held instruments (manipulators, forceps, biopsy punch).

The large-diameters blade supports (9.50 and 10.00mm) - create one circular incision in order to find cleavage plan using a manipulator. Forceps can then be used to peel the corneal endothelium.

Eventually, a biopsy punch can be used for the S-marking (#20034) technique to create a hole (centered on cardinal holes) in the exposed stroma.

The small-diameters blade supports (7.50, 7.75, 8.00, 8.25 and 8.50 mm) - allow for the completion of the corneal graft diameter, after having peeled the corneal endothelium, without being fully penetrating.

4-PRE-REQUISITES BEFORE USE AND DIRECTIONS FOR USE

To prepare device: follow step 1 to 4.

To maintain cornea with the base: follow step 5 to 10.

To proceed trephination with the blade support: follow step 11 to 12.

To finish the procedure: follow 13 to 14.

Preparing of device

- 1- Check the blade diameter which is marked on top of the blade support.
- 2- Remove the sterile blister from the pouch.
- 3- Remove the base, blade support and the syringe from the blister.
- 4- Check the integrity of the blade under a microscope.

Maintaining the cornea

- 5- Connect the tubing to the syringe with the Luer lock connection (see Fig. 3).
- 6- Press the piston of the syringe and hold.
- 7- Place the graft on the cutting block base, endothelial side up.
- 8- Center the graft using the groove centering ring
- 9- Release the piston of the syringe to create the vacuum and hold the graft in position on the base. Check that the graft is still well centered.
If the graft is not well centered at the end of this preparation, repeat the steps 6, 7, 8 and 9 to obtain a perfect centering.
- 10- Proceed to trephination steps or other manual manipulation of the donor cornea

Trephination

- 11- Place the blade support on the base, let it lower down by gravity to avoid setting it with a wrong angle (risk of blockage) (see Fig. 4 and Fig. 5) then press down gently.

Warning: Do not turn the punch during or after trephination.

- 12- Gently pull the blade support upward from the base. The graft should stay in place on the base as a result of the vacuum. Then completely remove the blade support. (see Fig. 6)

End of procedure

- 13- Disassemble the syringe from the tubing to release the vacuum. The donor button is now ready to be used on the patient.
- 14- At the end of the surgical procedure, discard the ONE Single-Use Guarded Punch into an appropriate container for destruction (according to appropriate regulations).

5-WARNING / CAUTION / PRECAUTION / ADVERSE EVENTS

- This device must be used only by experienced ophthalmic surgeons or eye banks technicians.
- These instructions explain how to use this medical device only; and do not describe the steps necessary to perform a corneal graft. Manual manipulation of the donor cornea with hand held instruments, and the surgical procedure, remain the responsibility of the eye banks technicians and/or the surgeon experienced in performing keratoplasty.
- Do not use if expiration date has passed.
- Do not use if storage conditions have failed (humidity traces).
- Do not use if the packaging or device is damaged (check the integrity of pouch and inspect the blade under a microscope before use).
- Do not have any physical contact with the cutting edge of the blade to avoid injuries.
- Do not use in combination with other version of guarded and classic punch (series 17207Dyyy, 17206Dyyy, 17208Dyyy, 17200Dyyy).
- Stop the use of the device in case of any malfunction during the procedure, including potential complication described in part 6.
- The design of this medical device is 'single use only': reuse of this device may affect its mechanical and biological features, and its clinical performances; and may cause device failures or may exposes the patient to adverse events such as contamination, bacterial infections or allergic reactions.

6-POSSIBLE SIDE EFFECTS AND INFORMATION TO RELAY


To date, the identified risks by the manufacturer associated to punches are:

- Loss of suction
- Graft attachment on the blade
- Incomplete cut

Any serious incident occurring in conjunction with this device must be reported to the manufacturer and the competent authority of the user's and/or patient's Member State.

7-LIST OF REFERENCES

Reference: #17212Dxxx	
17212D750	ONE-DEEP WELL GUARDED PUNCH D7.50 MM X1
17212D775	ONE-DEEP WELL GUARDED PUNCH D7.75 MM X1
17212D800	ONE-DEEP WELL GUARDED PUNCH D8.00 MM X1
17212D825	ONE-DEEP WELL GUARDED PUNCH D8.25 MM X1
17212D850	ONE-DEEP WELL GUARDED PUNCH D8.50 MM X1
17212D950	ONE-DEEP WELL GUARDED PUNCH D9.50 MM X1
17212D1000	ONE-DEEP WELL GUARDED PUNCH D10.00 MM X1
Reference: #17213Dyyy	
17213D950	ONE-DEEP WELL GUARDED COVER D9.50 MM X1
17213D1000	ONE-DEEP WELL GUARDED COVER D10.00 MM X1

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Moria Inc