EN - M2, M2SU

INSTRUCTION MANUAL 1-13

FR - M2, M2SU

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ES - M2, M2SU

MANUAL DE INSTRUCCIONES 29-41

IT - M2, M2SU

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DE - M2, M2SU

BENUTZERHANDBUCH..... 57-69

NL - M2, M2SU

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PT - M2, M2SU

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SV - M2, M2SU

INSTRUKTIONSMANUAL 99-111

TR - M2, M2SU

KULLANIM TALIMATLARI
EL KİTABI 113-125

Ref.: 65013

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MORIA S.A. – 15, rue GEORGES BESSE - 92160 ANTONY – France
#65013-F-02.2017

CONTENTS

1	DISCLAIMER	2
1.1	CARE AND MAINTENANCE	2
1.2	USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES	2
1.3	HEAD CALIBRATION VALUE AND FLAP THICKNESS	2
2	EQUIPMENT LIST	3
3	REGULATORY INFORMATION	4
4	LABELING INFORMATION	4
5	IMPORTANT INFORMATION	5
5.1	DESCRIPTION	5
5.2	INDICATIONS	5
5.3	CONTRAINDICATIONS	5
5.4	WARNINGS	5
5.5	POTENTIAL ADVERSE EVENTS	5
5.6	PRECAUTIONS	6
5.7	PRODUCT COMPLAINTS	6
5.8	PACKAGING	6
6	INSTALLATION AND CONNECTION	7
7	OPERATION	8
8	TROUBLESHOOTING	10
9	CARE & MAINTENANCE	10
9.1	INTRODUCTION	10
9.2	INITIAL TREATMENT AND STORAGE	10
9.3	CLEANING-DISINFECTION-STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE	10
9.4	CLEANING-DISINFECTION-STERILIZATION-STORAGE OF THE SUCTION RINGS AND M2 METALLIC HEADS	11
10	WARRANTY	12
10.1	SCOPE OF WARRANTY	12
10.2	NON-APPLICATION OF WARRANTY	12
10.3	WARRANTY PERIOD	12
10.4	LIABILITY	12
11	DRAWINGS	13

1 DISCLAIMER

1.1 CARE AND MAINTENANCE

The M2-M2SU system has been designed for optimal operation, provided that the recommendations listed in this user manual are followed carefully. If, for any reasons, the system does not perform properly, have it checked immediately by MORIA. In order to maintain the original performance of your microkeratome, MORIA strongly recommends annual maintenance and servicing of all its reusable products.

As only MORIA and its agents are fully expert in MORIA products, servicing and maintenance must be carried out by MORIA or its approved agents.

MORIA shall not be held liable for any malfunction or damage to the apparatus, poor results, or surgical complications due to maintenance being having been carried out by an unqualified operator or third party.

Any such unauthorized intervention shall render the guarantee and any maintenance contract null and void.

1.2 USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES

The materials used in the microkeratome blade, blade holder, and head have been selected for their biocompatibility and sliding characteristics. The dimensions and tolerances of the blade have been determined in consideration of the dimensions and tolerances of the head of the keratome. MORIA's manufacturing and inspection procedures guarantee there is no conflict in dimensions between head and blade, and that the blade will slide smoothly in the head.

Single-use devices should not be re-used. Doing so will negatively affect their clinical performance and increase the potential for adverse events.

The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and damage the microkeratome.

MORIA shall not be held liable in the event of a malfunction or damage to the microkeratome, poor results or surgical complications due to the reuse of a single-use product, or the use of consumables other than those supplied by MORIA.

MORIA handpieces must only be connected to MORIA devices (console unit, heads, suction rings, etc.).

All warranties become null and void if the microkeratome degrades or malfunctions due to such practices.

1.3 HEAD CALIBRATION VALUE AND FLAP THICKNESS

Flap thickness is a key factor for LASIK surgery. Multiple parameters affect flap thickness and standard deviation. Numerous scientific studies have demonstrated that flap thickness is influenced by several patient-related factors such as keratometry (K) readings, corneal anatomy, preoperative pachymetry and refraction, intraocular pressure (IOP), and also by surgery-related factors, such as corneal hydration and the speed of the pass (when using manual microkeratomers).

Ultrasonic pachymetry measurements are not always accurate and reproducible, and results may also vary with surgeons' techniques and calibration of the device. Thus, a head labeled 130 cuts on average a 160-micron flap:

- 130 microns is a dimension on the head that can be measured and which corresponds to a distance from the cutting edge of the blade to the front plate of the head
- 160 microns is the average flap thickness that will be cut by this head, with variations around this average value.

In the example above, the "incremental value" for a 130 head is 30 microns on average.

Head calibration value	Incremental value	Average flap thickness
90/110/130	+30	120/140/160

MORIA can only guarantee the dimensional characteristics of the head, not the surgical result.

For the EVOLUTION 2, 3 and 3E consoles, please refer to the related instruction manuals (#65016, 65038, 65051, 65060).

The most recent version of this user guide and additional information on your keratome are available on the MORIA website: <http://www.moria-surgical.com>.


2 EQUIPMENT LIST

Designation	MORIA reference	System	
		M2	M2SU
M2-M2SU suction ring (size -1, 0, +1, +2, +3)	19325/X X = -1/0/1/2/3	X	X
M2-M2SU suction ring size -1 Large-Cut	19379/-1	X	X
M2 metallic head	19327/110 19327/130	X	
Box of 10 sterile disposable M2 blades	19329	X	
Box of 10 sterile disposable M2SU heads	19334/90 19334/130		X
M2 motor	19326	X	X
Sterilisation box	22519330	X	X
Storage box	22519331	X	X
Tonometer	19042	X	X
Aspiration tubing	19138	X	X
Dual speed M2 connecting cable (2 forward speed)	19335	X	X
EVOLUTION 2 console	19350	X	X
EVOLUTION 3 console	19360	X	X
EVOLUTION 3E console	19380	X	X
EVOLUTION 2 footswitch	19351	X	X
EVOLUTION 3-3E footswitch	19361	X	X
Footswitch Epi-K™	19381	X	X
M2-M2SU user manual	65013	X	X
EVOLUTION 2 console user manual	65016	X	X
EVOLUTION 3 console user manual	65038	X	X
EVOLUTION 3E console user manual (serial number below 5000)	65051	X	X
EVOLUTION 3E console user manual (serial number 5000 and above)	65060	X	X
Annexe "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	65073	X	X
M2-M2SU nomogram (ES)	65019	X	X
M2-M2SU nomogram (FR)	65023	X	X
M2-M2SU nomogram (EN)	65024	X	X
M2-M2SU nomogram (XX=other)(MORIA website)	65024XX	X	X
M2-M2SU nomogram (USA only)	65031	X	X

3 REGULATORY INFORMATION

	MORIA S.A. 15, rue Georges Besse - 92160 Antony - France Phone +33 (0) 146 744 674 Fax +33 (0) 146 744 670 moria@moria-int.com http://www.moria-surgical.com
Customer Service Information	Contact your local dealer or MORIA
EUROPE	CE 0459 according to MDD 93/42/CEE
USA	Product registered at the Food and Drug Administration (FDA): M2SU K 022560 Caution USA only: Federal law restricts the use of this device to physicians or licensed practitioners.
Electrical Safety Standard	IEC 60601 - Classe II BF
	For EU customers only: this symbol indicates that within the European Union, the product must be discarded in a separate collection bin at the end of its useful life. This applies not only to this device, but also all accessories, including the footpedal and electrical motors, regardless of whether those accessories are marked with the symbol. Do not dispose as unsorted municipal waste. For users out of European Union: please refer to local environmental regulatory regarding waste of electrical and electronic equipment.

4 LABELING INFORMATION

	CATALOGUE REFERENCE
	QUANTITY
	USE BY
	BATCH CODE
	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
	MANUFACTURER
	DATE OF MANUFACTURE
	CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)
	CONSULT OPERATING INSTRUCTIONS
	DO NOT USE IF PACKAGE IS DAMAGED
	STORAGE CONDITIONS: <ul style="list-style-type: none"> • TEMPERATURE: XX°C – YY°C / XX°F – YY°F • HUMIDITY RATE: XX% – YY%
	NON-STERILE
	KEEP DRY
	DISCARD IN A SEPARATE COLLECTION BIN
	CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.

5 IMPORTANT INFORMATION

5.1 DESCRIPTION

The M2-M2SU microkeratome is an automated, mechanical, pivoting microkeratome. It has three components:

- a M2 motor (#19326) containing 2 independent motors: one for advancement (rotation), and one for oscillation
- a metallic suction ring (#19325/xx, 19379/-1)
- a single use plastic head (M2SU) encasing a pre-inserted blade, or a metallic reusable head (M2).

The M2 microkeratome operates with the EVOLUTION 2, 3, and 3E control units (#19350, 19360, 19380).

Please refer to the user guides:

- #65016 (EVOLUTION 2)
- #65038 (EVOLUTION 3)
- #65051 (EVOLUTION 3E serial numbers below 5000)
- #65060 (EVOLUTION 3E serial numbers 5000 and above).

5.2 INDICATIONS

The M2 microkeratome is indicated for the creation of lamellar corneal flaps on patients' eyes when using M2-M2SU heads in combination with a suction ring in corneas with preoperative pachymetry of 500 microns or greater, and keratometry between 39 D and 49 D.

During laser in-situ keratomileusis (LASIK), the flap is then lifted to enable photoablation of the stroma with an excimer laser.

5.3 CONTRAINDICATIONS

Patients who are not candidates for LASIK.

Note: Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to determine the risk/benefit ratio. Particular attention must be given before performing a keratectomy on a patient with any of the following conditions:

- preoperative pachymetry under 500 microns
- keratometry under 39 D
- keratometry over 49 D
- inability to withstand a transient rise in intraocular pressure.

5.4 WARNINGS

- **Single-use devices should not be re-used. Doing so will negatively affect their clinical performance and increase the potential for adverse events.**
- Do not use disposable materials and/or components of a brand other than MORIA with the M2 microkeratome.
- Blades must only be inserted in heads by hand. Never use any tools or spanners. Incorrect assembly may cause incomplete or uneven cuts due to lack of blade oscillation.
- The heads must only be screwed on by hand. Never use any tools or other spanners. Incorrect assembly may cause incomplete or uneven cuts due to lack of blade oscillation.
- Never pull on the cable connected to the handpiece, and never hold the motor by the cable.

5.5 POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. LASIK surgery requiring the use of a microkeratome that cuts a corneal flap, potential side effects of laser refractive surgery may include but are not limited to: visual anomalies, dry eye and flap related complications (free cap, incomplete flap, buttonhole, epithelial defect, flap dislocation, flap striae, wrinkles, etc.).

Inappropriate use, deterioration of microkeratome and/or non-respect to contraindications (§ 5.3) and warnings (§ 5.4) expose patient at higher risk to adverse events.

Note: Additional surgery may be necessary to correct some of these potential adverse events.

5.6 PRECAUTIONS

- The keratotomy must be performed only by experienced refractive surgeons with specific training in the use of the M2-M2SU microkeratome.
- Preoperative and operative procedures, including knowledge of surgical techniques, proper head and ring selection, and assembly and placement of the microkeratome are important considerations in the successful use of the system by the surgeon. Furthermore, the proper patient selection and compliance will greatly affect the results.

For USA only: CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

Other preoperative, intraoperative, and postoperative warnings and precautions are as follows:

• PREOPERATIVE

- o Only patients that meet the criteria described in the indications should be selected.
- o Care should be used in the handling and storage of the microkeratome components. They should not be scratched or otherwise damaged. The M2 motor should be protected during storage, especially from corrosive environments.
- o Check the label and expiration date on unopened packages of M2SU heads and M2 blades.
- o After opening the sterile package of M2SU heads, verify that head information is consistent with information on outer package labeling.
- o The head, suction ring, M2 motor, and the control unit should be fully inspected prior to use. The blade should be inspected on both sides under a microscope.
- o The surgeon should be familiar with the various components before using the microkeratome and verify that all parts and necessary instruments are present and properly assembled before the surgery begins. Additional sterile components should be available in case of an unexpected need.
- o The selection of the proper head, ring, and stop position for each eye is crucial to the success of the procedure: see MORIA nomogram.

• INTRAOPERATIVE

- o Breakage, slippage, or misuse of the microkeratome and its components may cause injury to the patient or operating personnel.
- o Before any usage, lubricate the eye, the suction ring, the head, and the blade with a physiological saline solution or another appropriate ophthalmic solution that is compatible with the M2-M2SU components.

• POSTOPERATIVE

- o The surgeon's postoperative instructions to the patient and the corresponding patient compliance are extremely important.

5.7 PRODUCT COMPLAINTS

Health care professionals (e.g., customers or users of this device) having any complaints about or dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify MORIA or its distributor by telephone, fax or written correspondence, and have the products checked by MORIA.

When filing a complaint, please provide the component(s) name(s), reference(s), batch number(s), as well as your name and address, the nature of the complaint, and the relevant patient data. Also disinfect and return the components.

5.8 PACKAGING

Packages of each of the components should be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness, and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used and should be returned to MORIA.

Do not use single-use products if package is open or torn.

6 INSTALLATION AND CONNECTION

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Select the head	<ul style="list-style-type: none"> • M2: select the appropriate pre-calibrated head. The microkeratome head is labeled on its front side and indicates depth of cut. • M2SU: select the appropriate pre-calibrated head. The microkeratome head is labeled on its front side and indicates depth of cut. Disposable M2SU heads are supplied sterile and for single use only. Check that package is undamaged, unopened and within the expiration date. 	A, B
2a	M2: Insert the blade into the blade housing of the head	<ul style="list-style-type: none"> • Lubricate the M2 head with balanced salt solution or another appropriate ophthalmic solution. • Inspect the blade and particularly the cutting edge before and after insertion. Blades with spots, nicks, debris or other imperfections which could compromise performance should be discarded. • Insert the blade into the blade housing of the head, being careful not to touch the front edge of the blade. • Press gently on the blade and slide it completely into the head. In case the blade does not slide smoothly into the head, check for debris or deposits in the blade housing or on the blade itself. If, after careful checking of the blade housing, it does not slide, replace the blade with a new one. • Guide the blade into the head avoiding contact of the cutting edge with metal part of the head. • The blade design allows for only one direction for insertion and prevents incorrect assembly. • Before mounting the M2 head onto the turbine motor, make sure that the blade holder slot is centered within the head (centered within the female collar). 	C, D
2b	M2SU: lubricate the head.	<ul style="list-style-type: none"> • M2SU heads are supplied sterile with pre-inserted blade to avoid blade manipulation. • Inspect the blade and particularly the cutting edge before and after insertion. Blades with spots, nicks, debris or other imperfections which could compromise performance should be discarded. • Place the protective case on a flat surface and lift the top of package with caution. • Fill the cavity of the plastic package with balanced salt solution or another appropriate ophthalmic solution to fully wet the head. 	B, E
3	Assemble the head into the M2 motor	<ul style="list-style-type: none"> • Before using the M2 electrical motor ensure that the motor is clean and in proper working order. • Two different M2 motors exist: with and without the locking lever. • M2: the electrical motor locks onto the keratome head: <ul style="list-style-type: none"> o holding the microkeratome head, align the flat areas of the head and the motor o turn the head 1/4 turn counter-clockwise onto the nose of the motor (finger tight only). o check that the head is securely locked and tightened to the motor • M2SU: <ul style="list-style-type: none"> o keep the head in the package while setting it onto the motor. o as for the reusable M2 head, screw the nose of the electrical motor approximately 1/4 turn counter-clockwise onto the head o keep the blister aside • Check that the M2-M2SU head is securely locked to the electrical motor. • If the head is not completely screwed on, the electrical motor will not drive the blade holder and the blade will not move. Poor or irregular resection can result from incorrect assembly. 	F, G, H, I
4	Inspect the head after assembly	<ul style="list-style-type: none"> • The head must be perfectly clean inside and outside and free from debris, particles, oxidation and deposits. There should be no scratches on the plates. If any, replace the head. • After assembly and using a microscope, carefully inspect the blade under high magnification to ensure that it is not damaged. If so: discard it (M2 blade, M2SU head) in an appropriate container and replace it. 	

5	Check for blade oscillation	<ul style="list-style-type: none"> Refer to EVOLUTION 2, 3 or 3E user manual (#65016, 65038, 65051, 65060) for operations prior use. Activate M2 motor to check for smooth, strong and uninterrupted blade oscillation. Do not use the M2 motor if the blade oscillation is not smooth, regular and uninterrupted. 	
6	Select suction ring according to nomogram, then inspect it	<ul style="list-style-type: none"> Refer to the indicative nomogram for the selection of the speed, ring size and stop value. The suction ring must be carefully inspected prior to use; it must be perfectly clean, free from debris, particles and deposits. Carefully check that the suction port is not obstructed. The hinge will be placed facing the arrow on the suction ring. 	J
7	Set the stopper	<ul style="list-style-type: none"> The adjustable stop pin is pre-mounted on the suction ring and used to determine the size of the hinge. Two types of stopper exist depending suction ring: « 3 stop position» or « 4 stop position». To configure the stop, follow these steps: <ul style="list-style-type: none"> unscrew the stop setting screw lift it up then turn the stop to the appropriate value the value of the stop selected must line with the etched ring size screw the stop setting screw. 	K, L
8	Connect the suction ring to the aspiration tubing	<ul style="list-style-type: none"> The aspiration tubing (#19138) is supplied sterile and for single use only. Check that package is undamaged, unopened and within the expiration date. Inspect the aspiration tubing and in case of kinks or obstructions, replace it. 	M

7 OPERATION

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Preliminary check before any usage	<ul style="list-style-type: none"> It is imperative to verify that the devices are perfectly clean and free of any organic or other residues. 	
2	Connect the M2 motor to the control unit. Check vacuum prior use	<ul style="list-style-type: none"> See EVOLUTION 2, 3 or 3E user manual (#65016, 65038, 65051, 65060) for operations prior use. 	

3	Assemble the M2 motor in the suction ring	<ul style="list-style-type: none"> • Before locking the assembled M2 microkeratome motor and M2 head onto the suction ring, make sure that the rotation shaft of the motor is in contact with the rear stop. • If the rotation shaft is not in the proper position, the M2 will not be oriented correctly on the suction ring and will not properly lock. • If the rotation shaft is not in the correct position, activate the vacuum and then the backward (reverse direction) footswitch to correctly position the rotation shaft. • The M2-M2SU head pivot opening should fit easily over the pivot post of the M2 suction ring. • Align suction ring pivot post and M2 head pivot hole. Carefully place the assembled M2 keratome and M2 motor onto the suction ring. • Avoid inadvertent contact with suction ring pivot post and M2 head blade area. If contact occurs, blade edge may become damaged, and a new blade (M2 head) or head (M2SU head) should be used. • Engage the microkeratome head on the suction ring pin. • Turn the locking ring clockwise until it is locked: <ul style="list-style-type: none"> o ring should lock smoothly without requiring force o two black marks on the ring become hidden and the third is aligned with a black mark on the motor. • WARNING: DO NOT START THE CUT IF THE LOCKING RING IS NOT CORRECTLY LOCKED. 	N, O, P
4	Check for blade oscillation	<ul style="list-style-type: none"> • Refer to EVOLUTION 2, 3 or 3E user manual (#65016, 65038, 65051, 65060) for operations prior use. • Activate M2 motor to check for smooth, strong and uninterrupted blade oscillation. • Do not use the M2 motor if the blade oscillation is not smooth, regular and uninterrupted. 	
5	Position the suction ring on the eye. Activate vacuum by pressing the "Vacuum" footswitch once.	<ul style="list-style-type: none"> • Verify that the device will not interfere with any other medical devices used during the procedure. 	
6	Check the intraocular pressure with the tonometer	<ul style="list-style-type: none"> • The tonometer (#19042) must be perfectly dry and used only on dry eyes. • If the pressure is below 65 mm Hg, do not proceed with the surgery. 	
7	Lubricate the suction ring, the head and the blade	<ul style="list-style-type: none"> • Prior to use, lubricate the rails of suction ring, head, and blade with balanced salt solution or another appropriate ophthalmic solution. • The device should not be placed in contact with solutions for which the potential interactions are unknown. 	
8	Once the system is in position, activate blade oscillation by pressing the "Forward" footswitch. When the microkeratome head touches the stop, immediately release the "Forward" footswitch. Reverse by pressing the "Backward" footswitch.	<ul style="list-style-type: none"> • See EVOLUTION 2, 3 or 3E console user manual (#65016, 65038, 65051, 65060). • Hold the suction ring by its handle; check that your fingers do not hinder the forward movement of the device. • Check nothing prevents or modifies head movement in the rails of the suction ring • Ensure that there are no obstructions (speculum, eyelashes, eyelids, etc.) in its path 	
9	Release the "Vacuum" footswitch by pressing once. Remove the suction ring.	<ul style="list-style-type: none"> • Discard the disposable blade (M2 head) or disposable head (M2SU head) and aspiration tubing after each procedure in an appropriate container. 	

8 TROUBLESHOOTING

For any further information and description, please refer to EVOLUTION 2, 3 or 3E console user manual (#65016, 65038, 65051, 65060).

9 CARE & MAINTENANCE

The recommendations for cleaning, disinfection and sterilization are updated and available on MORIA's website: <http://www.moria-surgical.com>.

9.1 INTRODUCTION

In case of unusual vibrations or noises, do not use the unit and contact your distributor.

Please contact MORIA for any other servicing.

The following are only recommendations. They must be adapted according to the laws in the country of utilization. All the products and solutions for cleaning and disinfection must be used according to the manufacturer's instructions.

MORIA recommends that you:

- avoid touching potentially contaminated areas,
- use gloves when performing cleaning and decontamination tasks.

9.2 INITIAL TREATMENT AND STORAGE

To avoid risks of condensation inside the packaging, unpack and store the devices in a clean, dry environment.

Do not store the instruments in an environment or next to other products that may possibly have a corrosive or magnetic effect. Avoid any contact between devices, especially those comprised of different materials.

The device should be damage free and have no scratches or other surface defects.

Fragile devices or those requiring particular manipulation must be handled separately, taking special care to protect delicate parts.

9.3 CLEANING-DISINFECTION- STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Cleaning <ul style="list-style-type: none">• Dip the blade oscillation shaft of the motor in a sterile distilled water bath and run the motor back and forth 10 times.• Remove the motor from the cup, unplug it from the console, and wipe the oscillation shaft of the motor with a lint-free cloth.• Cable can be cleaned with a lint-free cloth moistened with cleaning solutions.• Air-dry the shaft of the motor with clean, filtered compressed air (medical compressed air).	<ul style="list-style-type: none">• The handpiece must be carefully cleaned after each surgical procedure.• Do not use abrasives or scrapers to clean any elements of the system. Doing so may result in reduced precision, burrs and/or irregular tissue resection.• Do not immerse the motor in the sterile distilled water bath.• To prevent damaging the connectors and to ensure proper functioning of the motor, never pull on the cables and never hold the motor by its cable.	Q
2	Disinfecting & Drying <ul style="list-style-type: none">• Wipe each part of the motor with a lint-free cloth moistened with disinfecting solutions.• Use appropriate disinfecting solutions (spray or wipe clean disinfectant) according to the manufacturers' instructions.• Dry carefully with disposable instrument cleaner (lint-free) then with clean, filtered, compressed air.	<ul style="list-style-type: none">• It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.• Do not gas sterilize (ETO) the motor.• Do not autoclave the motor.• No current sterilization process is compatible with MORIA motors	
3	Storage	<ul style="list-style-type: none">• Failure to completely dry the inside of the motor can result in oxidation.• Do not store products that are not completely dry; doing so may create rust and result in irregular corneal resection.• When not in use, the M2 system should be kept in its storage box (#22519331) in a dry atmosphere.	

9.4 CLEANING-DISINFECTION-STERILIZATION-STORAGE OF THE SUCTION RINGS AND M2 METALLIC HEADS

Steps	What to do	 IMPORTANT WARNINGS
1	<p>Double cleaning</p> <ul style="list-style-type: none"> On leaving the operating theatre, immediately immerse the devices in an aldehyde-free detergent/disinfectant solution (e.g.: Alkazyme® / Alkapharm Laboratories) for at least 15 minutes according to manufacturer recommendations. Manual cleaning (or use a machine, provided that it does not recycle the cleaning products, after disassembling the device where appropriate). Rinse the circuit with water. Immerse the device in an aldehyde-free detergent/disinfectant solution (a different bath from the one above). Manual cleaning (or use a machine, provided that it does not recycle the cleaning products). Rinse using filtered distilled or osmosis-treated water. Wipe the device using a non-woven, lint-free, disposable cloth. 	<ul style="list-style-type: none"> Manual cleaning: mechanical cleaning of the medical device using soft brushes (metal brushes and scouring pads are prohibited). IMPORTANT: Change the baths after each use. Incompatibility: <ul style="list-style-type: none"> Do not use bleach Do not apply this protocol to products made of aluminium alloy. Do not apply this protocol to heat-sensitive (motor and turbine).
2	<p>Deactivation of NCTAs (Non-Conventional Transmitted Agents)</p> <ul style="list-style-type: none"> On leaving the operating theatre, immediately immerse the devices in an aldehyde-free detergent/disinfectant solution, preferably enzymatic (e.g.: Alkazyme® / Alkapharm Laboratories) for at least 15 minutes according to manufacturer recommendations, and then rinse the devices. Cleaning, disinfectant and rinsing fluids must be stored in sealed containers and disposed of by appropriate procedures for contaminated biological liquid waste. Immerse the device in a soda solution (1 N) for one hour. Rinse the device manually three times (using water from the network), and check the pH when performing the last rinse to prevent any risk of burns (neutral pH). The last rinse must be carried out using sterile filtered distilled or osmosis-treated water. Wipe the device using a non-woven, lint-free, disposable cloth, and reinstall the device in the normal disinfection/sterilisation circuit. 	<ul style="list-style-type: none"> Deactivation of NCTAs: for patients with a risk of Creutzfeldt-Jakob disease only – CJD questionnaire (see the French Ministry's guide, dated November 2005) Note: If the pH is not neutral, perform further rinses until an acceptable pH is reached (close to pH7).
3	<p>Sterilisation</p> <p>Sterilise the sterilisation tray using an autoclave (moist heat steriliser) at a temperature of at least 134°C (273.2°F) for 18 minutes.</p>	

Bibliography:

- March 2001 – Circular No. DGS/5C/DHOS/E2/138 dated 14 March 2001 on risk management in health establishments
- November 2005 - Guide from the Ministry on the treatment of medical devices for ophthalmology and contactology - <http://www.sante.gouv.fr>.

10 WARRANTY

10.1 SCOPE OF WARRANTY

Designation	MORIA references
M2-M2SU suction ring (size -1, 0, +1, +2, +3)	19325/-1 19325/0 19325/1 19325/2 19325/3
M2-M2SU suction ring size -1 Large-Cut	19379/-1
M2 metallic head	19327/110 19327/130
M2 motor	19326
Sterilisation box	22519330
Storage box	22519331
Dual speed M2 connecting cable (2 forward speed)	19335
EVOLUTION 2 console (except the battery)	19350
EVOLUTION 3 console (except the battery)	19360
EVOLUTION 3E console (except the battery)	19380
EVOLUTION 2 footswitch	19351
EVOLUTION 3-3E footswitch	19361
Footswitch Epi-K™	19381

- The above items as well as spare parts and labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, after having previously been disinfected.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.
- Any defective part exchanged during the period of warranty becomes property of MORIA.
- Instruments and accessories that cannot be reused are excluded from this warranty.
- Future upgrades and/or improvements on the keratome are not implied by this warranty.

10.2 NON-APPLICATION OF WARRANTY

The warranty will not be applicable under any of the following conditions:

- Defects or malfunction that occur out of the warranty period (10.3).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- The use of supplies, spare parts, or accessories other than those supplied by MORIA.
- Any disassembly, modification or intervention carried out on the devices by a person not authorized by MORIA.

10.3 WARRANTY PERIOD

- The warranty takes effect on the day the equipment is shipped.
- The warranty is for a period of 12 months from the effective date.

10.4 LIABILITY

- The liability of MORIA is limited to the supply of the services mentioned in paragraph 10.1. MORIA will not be held responsible for any direct or indirect damage suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre (France) will have sole jurisdiction.

11 DRAWINGS

A. M2 HEAD

- 1 Head calibration value
- 2 Female collar
- 3 Blade housing
- 4 Flap housing

B. M2SU HEAD

- 1 Blade

C. M2 BLADE

- 1 Cutting edge
- 2 Blade holder

D. M2 HEAD WITH INSERTED M2 BLADE

- 1 M2 head
- 2 M2 blade

E. OPEN PEEL POUCH

- 1 Place the protective case on a plane surface
- 2 Lift the top of the package with caution

F. LOCKING LEVER

- 1 M2 motor without locking lever
- 2 M2 motor with locking lever

G. DISTAL SHAFT OF M2 MOTOR (MOTION AND OSCILLATION)

- 1 Oscillation of the blade
- 2 Rotation of the keratome on the eye

H. MOUNTING THE M2 HEAD ON THE M2 MOTOR

- 1 Engagement of motor inside the head: align the 2 black dots
- 2 Locking of the head onto the motor: screw counter-clockwise the nose of the electrical motor onto the head

I. MOUNTING THE M2SU HEAD ON THE M2 MOTOR

- 1 M2 motor
- 2 M2SU head
- 3 Package of M2SU head

J. M2 SUCTION RING

- 1 Handle
- 2 Guide
- 3 Aperture
- 4 Adjustable stop ring
- 5 Pre-set starting and ending position mark
- 6 Identification of suction ring
- 7 Pivot
- 8 Mark of hinge position

K. CLOSE VIEW OF THE STOP RING ON M2 SUCTION RING

- 1 Safety screw
- 2 Stop ring

L. VIEW OF THE STOP RING

- 1 3 stop positions
- 2 4 stop positions

M. CONNECTING THE SUCTION RING TO THE ASPIRATION TUBING

- 1 Aspiration tubing end
- 2 Suction ring end

N. LOCKING THE M2 MICROKERATOME MOTOR AND HEAD ONTO THE SUCTION RING

- 1 Rotation shaft NOT in contact with rear stop
- 2 Rotation shaft in contact with rear stop

O. ENGAGEMENT OF HEAD ON THE M2 RING

- 1 M2 motor with M2-M2SU head
- 2 M2 suction ring

P. M2 LOCKING RING SYSTEM

Unlocked

- 1 Turn the locking ring until the 2 marks are aligned
- 2 Turn the locking ring until the 2 marks are aligned

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