# #65093

EN

## Instructions for use of the ONE® Muraine Vacuum Punch #17209Dxxx

### Indication for use

The Moria Muraine Vacuum Punch is designed to perform lamellar keratoplasty, in which 2 corneal incisions [I] (see Fig.1) are made at a depth of approximately 300µm.

17209Dxxx (Dxxx = blade diameter in mm)

### Device components

- Cutting block base (A) (see Fig.2)
- Trephine with a blade of which the descent is limited (B) (see Fig.2)
- Syringe (C) (see Fig.2)

### Warning

This instruction for use only refers to the use of medical device but does not describe surgical steps to perform keratoplasty. Graft surgery must be performed only by experienced corneal surgeons.

- Check the integrity of the pouch.
- Check the expiration date on the labels.
- Inspect the blade under a microscope.
- The device is sterile and ready for use. Do not use it if the pouch in which it is packaged is damaged.
- **The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and damage the device.**
- **Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performances and exposes the patient to adverse events.**

Instructions for use and any updates are available on Moria (<http://www.moria-surgical.com>) website.

The most recent version of this document is available at <http://www.moria-surgical.com> via the restricted access heading «Documentation/Library». Translations are also available in other languages under that same heading. Registration required.

### Directions

- 1 / Check the blade diameter stamped on top of the trephine blade housing.
- 2 / Remove the sterile blister from the pouch.
- 3 / Remove the cutting block base, trephine blade housing, and the syringe from the blister.
- 4 / Connect the tubing to the syringe with the Luer-Lock™.
- 5 / Check the integrity of the trephine blade under the microscope (see Fig.3).
- 6 / Press the piston of the syringe half way and hold.
- 7 / Place the tissue on the cutting block base, endothelial side up. Center the graft using the 4 cardinal holes
- 8 / Release the piston of the syringe to create a vacuum that will hold the graft in position on the cutting block base.  
Check that the graft remains well centered.  
If the graft is not well centered at the end of this preparation, repeat steps 6 to 8 to obtain precise centering.
- 9 / Place the trephine blade housing on the cutting block base and press down firmly to cut the graft.  
**Warning: Do not turn the trephine during or after trephination.** (see Fig.4)
- 10 / Gently pull the trephine housing blade up from the cutting block base.  
The tissue should be kept in place on the cutting block base due to aspiration. Then completely remove the trephine housing blade.
- 11 / Disconnect the syringe to release vacuum.  
The tissue can now be placed on the Single-Use Artificial Chamber (Moria #19182) endothelium side up to complete the delamination.
- 12 / At the end of the surgical procedure, discard the punch in an appropriate container for destruction (according to the regulations in force).

### Accessories

ONE® – Artificial Chamber #19182 – Packaging : 1 per box, sterile, single-use.

### Patent & licences


US : Patent Pending

This device is manufactured and sold worldwide by Moria under an exclusive license from Pr. M. Muraine and University Hospital of Rouen, France.

**CAUTION for USA only: US federal laws restrict the sale of this device to or on the order of a physician.**

Explanation of symbols (General list of symbols used by Moria)

|                                                                                              |                                                                                                                                                 |
|----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
|  XXXXXX     | CATALOGUE REFERENCE                                                                                                                             |
|             | QUANTITY                                                                                                                                        |
|  YYYY/MM/DD | USE BY                                                                                                                                          |
|  XXXXXX     | BATCH CODE                                                                                                                                      |
|             | STERILIZED USING ETHYLENE OXIDE                                                                                                                 |
|             | DO NOT REUSE                                                                                                                                    |
|  XXXXXXX    | MANUFACTURER                                                                                                                                    |
|  YYYY/MM/DD | DATE OF MANUFACTURE                                                                                                                             |
|  XXXXXX     | CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)                                                                                                      |
|  XXXXXX     | CONSULT OPERATING INSTRUCTIONS                                                                                                                  |
|             | DO NOT USE IF PACKAGE IS DAMAGED                                                                                                                |
|             | STORAGE CONDITIONS: <ul style="list-style-type: none"><li>• TEMPERATURE: 15°C – 35°C / 59°F – 95°F</li><li>• HUMIDITY RATE: 30% – 70%</li></ul> |
|             | KEEP DRY                                                                                                                                        |
|             | DISCARD IN A SEPARATE COLLECTION BIN                                                                                                            |
|             | CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.                                         |

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**Moria Inc**
(CHINA) CO., LTD.**