
#19175

EN

Addendum for user manual #65009 & Instruction for use #65090

Moria informs the availability of the One Use-Plus (OUP) handpiece for artificial chamber (#19175) used with One Use Large-Cut heads (OU LC) (#19184/xxx) for single use or reusable artificial chamber.

1- Indications

OU LC heads for artificial chamber (#19184/xxx) guided for:

- RU (reusable) artificial chamber by the OUP guide ring for ALTK (#19173) mounted on ALTK artificial chamber (#19161) equipped with its universal cover (#19162)
- SU (single-use) artificial chamber by the OUP LC single-use guide ring (#19187) mounted on single-use artificial chamber (#19182)

enable to create, with an automated motion, a complete lamellar resection (without hinge) for donor graft preparation during keratoplasty surgery.

2- Reference & Designation

- **Reference:** #19175
- **Designation:** OUP handpiece for artificial chamber

3- Instructions for use

- 1/ Prepare Evolution 3E control unit (see instruction for use #65090 included in box of OU LC heads, and user manual Evolution 3E #65060/INTL):
 - perform vacuum test on Evolution 3E control unit
 - select OUP mode and speed 2
 - connect OUP handpiece to Evolution 3E control unit (Fig. 1)
- 2/ See instruction for use #65090 for preparation of artificial chamber, setting of the graft and preparation of the head.
- 3/ OUP handpiece setting:
 - screw the head onto the handpiece (Fig. 2)
 - inspect the blade under microscope and start the motor to check proper blade oscillation
 - set OUP handpiece on OUP guide ring #19173 or #19187 depending on artificial chamber used (Fig. 3)
- 4/ See instruction for use #65090 for cutting of the graft and end of procedure.

4- Maintenance

For the maintenance of the OUP handpiece for artificial chamber (#19175), please refer to user manuals #65039 or #65040 (§9.3).

As for the cleaning-disinfection-sterilization-storage of OUP handpiece: see chart on following page.

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Cleaning - To clean the screw core of the motor and metallic parts, use the disposable brush (#19149) moistened with cleaning solutions. - Dip the blade oscillation shaft of the motor in a sterile distilled water bath and run the motor back and forth 10 times. - Remove the motor from the cup, unplug it from the console, and wipe the oscillation shaft of the motor with a lint-free cloth. - Cable can be cleaned with a lint-free cloth moistened with cleaning solutions. - Air-dry the shaft of the motor with clean, filtered compressed air (medical compressed air).	The handpiece must be carefully cleaned after each surgical procedure. - Do not use abrasives or scrapers to clean any elements of the system. Doing so may result in reduced precision, burrs and/or irregular tissue resection. - Do not immerse the motor in the sterile distilled water bath. - To prevent damaging the connectors and to ensure proper functioning of the motor, never pull on the cables and never hold the motor by its cable.	Fig.4
2	Disinfecting & Drying - Wipe each part of the motor with a lint-free cloth moistened with disinfecting solutions. - Use appropriate disinfecting solutions (spray or wipe clean disinfectant) according to the manufacturers' instructions. - Dry carefully with disposable instrument cleaner (lint-free) then with clean, filtered, compressed air.	It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.	
3	Sterilization (optional) The OUP handpiece for artificial chamber is compatible with low temperature sterilization using hydrogen peroxide ("plasma sterilization"). Refer to the instructions for use provided by manufacturers of plasma sterilization.	Do not sterilize the motor with ethylene oxide (EtO). Do not autoclave the motor.	
4	Storage	- Failure to completely dry the inside of the motor can result in oxidation. - Do not store products that are not completely dry; doing so may create rust and result in irregular corneal resection. - When not in use, the OUP system should be kept in its storage box (#22519174) in a dry atmosphere.	

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MORIA Inc