#65113

EN

Instruction for use of ACP tubing set (#19192) ONLY FOR USE WITH ACP SYSTEM

IMPORTANT: Check the integrity of the pouches, as well as the expiration date on the labels. The medical device is sterilized with ethylene oxide and is for single use only. **THE DESIGN OF THESE MEDICAL DEVICES DO NOT ALLOW REUSE.** They should not be resterilized with ethylene oxide nor reused.

Contents


ACP tubing set can be divided in 2 parts:

- ① Air delivering tubing: one connector, one short tubing (5.12 in - 130 mm), one sterilizing hydrophobic filter (0.2 μ m), one long tubing (86.61 in - 2200 mm).
- ② Infusion tubing: one 3-way stopcock, one 20 mL syringe, one amortizing line (11.02 in - 280 mm), one clamp.

RECOMMENDATIONS:


- Before any connection, ACP system must be switched on (refer to user manual #65112).
- Remove all the sterile parts from the pouches. Place the base of the artificial chamber on the table and make sure it is stable.
- Use corneas with a rim having a minimum diameter of 16 mm and a well-centered cornea.

Use of ACP tubing set with an artificial chamber


Connect ACP tubing set:


- Connector: to ACP system output.
- Other extremity: to the artificial chamber.


Fill the syringe with corneal storage medium or a physiologic saline solution (BSS® or similar). Attach the syringe to the 3-way stopcock. Visually check that the 3-way stopcock is closed to ACP system, and open to both the syringe and the artificial chamber.

Infuse while positioning the cornea on the artificial chamber. If desired, apply cohesive viscoelastic solution to the endothelium. If needed, slightly lift the cornea with a pair of forceps to remove the possible air bubbles.

Precisely center the cornea on the artificial chamber and place the cover carefully in position. Check that the cornea is well centered. Seal the artificial chamber.


"OPEN system" = pressure directly coming from ACP system to the artificial chamber.

Before operation:

- Press on the footswitch to open the output valve (green LED ON - yellow LED OFF) and to release air in ACP tubing set (1).
- Shift the stopcock of ACP tubing set to be open to both ACP system and artificial chamber, while being closed to the syringe (2).
- Pressure is established in the artificial chamber. Corneal tissue is ready for operation.
Note: The pressure displayed on ACP system is the same as established in the artificial chamber.
If desired, check the pressure with a tonometer.

After operation:

- Shift the stopcock of ACP tubing set to be open to both artificial chamber and syringe, while being closed to ACP system (3).
- Press on the footswitch to close the output valve (yellow LED ON - green LED OFF) (4).

"CLOSED system" = pressure previously established with ACP system.

Before operation:

- Press on the footswitch to open the output valve (green LED ON - yellow LED OFF) and to release air in ACP tubing set (1).
- Shift the stopcock of ACP tubing set to be open to both ACP system and artificial chamber, while being closed to the syringe (2). Pressure is established in the artificial chamber.
Note: The pressure displayed on ACP system is the same as established in the artificial chamber.
- Shift the stopcock of ACP tubing set to be open to both ACP system and syringe, while being closed to the artificial chamber (3).
- Corneal tissue is ready for operation.
Note : Once the system is closed, syringe can gently be pushed away.
If desired, check the pressure with a tonometer.

After operation:

- Shift the stopcock of ACP tubing set to be open to both artificial chamber and syringe, while being closed to ACP system (4).
- Press on the footswitch to close the output valve (yellow LED ON - green LED OFF) (5).

Once operation is completed, unscrew the artificial chamber.

To prevent a possible collapse of the cornea, open the system to have sufficient counter-pressure, and remove the cornea.

At the end of the procedure, discard ACP tubing set into an appropriate container for disposal (according to the regulations in force).

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MORIA Inc