


Keratoplasty



*The right instruments whatever your
DMEK preferred technique*



**DMEK with
Moria**

**Preparing the donor graft
Preparing the recipient eye**




DMEK in 2021: a myriad of techniques

DMEK (Descemet Membrane Endothelial Keratoplasty) is a posterior lamellar keratoplasty technique that involves replacing the patient's damaged endothelium. As a minimally invasive technique, DMEK offers clinical benefits including rapid visual recovery¹⁻³ and low occurrence of rejection⁴. It has **therefore become a reference technique in endothelial keratoplasty**.


However, success relies on very delicate handling of the graft **to preserve as much as possible the endothelium and guarantee its viability**. Objective is to isolate the Descemet membrane with the endothelium to obtain a (purely) endothelial graft (without any posterior stroma). The difficulty of this technique led operators, eye banks and surgeons to innovate and standardize this procedure. That's why DMEK graft preparation technique has undergone **several evolutions**.


Reproducible DMEK grafts can now be obtained by choosing the most suitable harvesting technique for each specific user and setting among these 3:

- **manual dissection technique**
- **pneumatic dissection technique**
- **hydrodissection technique**



At Moria, we understand that every user has his preferred technique, so **we have developed a wide range of reusable and single-use instruments and devices to enable DMEK grafts to be performed using any technique** or, indeed, variation of a technique. In addition to being comprehensive, Moria range is also notable for its extremely high quality. Our reusable instruments in particular are **renowned worldwide for their durability and resistance**.






Preparing the donor graft


Common instruments whatever the harvesting technique

Holding forceps are used to handle a DMEK graft. Moria developed a range of fine toothed forceps dedicated to such a delicate step.




7835

Bonn forceps
0.1mm Micro-teeth
4mm Platforms
Reusable




7850A

Bonn-Moria forceps
0.1mm Micro-teeth
5mm Platforms
Reusable




13161

Bonn-Moria forceps
0.1mm Micro-teeth
3mm Platforms
Reusable




13160

Bonn forceps
0.1mm Micro-teeth
4mm Platforms
Reusable



17504X5
17504X10

Bonn forceps
0.12mm Micro-teeth
5mm Platforms
Single-use ☒



17221X10

Bonn forceps
0.12mm Micro-teeth
5mm Platforms
Single-use ☒



M1809

Bonn forceps
0.12mm Micro-teeth
Reusable




Preparing the donor graft

1) Manual dissection

The manual dissection technique for preparing DMEK grafts was introduced by Dr. Melles and his team in 2006¹ and has since evolved considerably. This technique consists in making a peripheral dissection at the trabecular meshwork, separating the Descemet membrane from the stroma, then peeling the membrane in balanced saline solution. Developments include the "Submerged Cornea Using Backgrounds Away" (**SCUBA**) technique⁵, the "**no-touch**" technique⁶, with use of a **guarded punch**^{7,8}, and the **SHARP** technique.^{9,10}

Moria offers a range of instruments and devices that will enable you to perform these various manual dissection techniques.




Instruments for peripheral dissection and cleavage plane identification

Peripheral dissection can be carried out using an angled hook, as is the case with the SHARP technique^{9, 10}, with a Y hook and micro-dissector as Dr. Price does¹¹, or using other types of dissectors.



20022

Micro-dissector
Active part 3mm
Reusable



20021

Hook
Descemet stroma hook
Reusable



6062A

Strampelli knife
Active part 3mm
Reusable

Preparing the donor graft

Endothelial peeling forceps

Peeling the endothelium requires one or two forceps² with flat, non-serrated jaws¹².




20038

Curved DMEK forceps

7° curvature, 10 mm platforms,
12 cm long.

Reusable



20039

Straight DMEK forceps

10 mm platforms,
12 cm long.

Reusable



17521x5
17521x10

Graft preparation forceps

1.2 mm oval-shaped tip
Curved, 15 mm long
Sold in box of 5 or 10 units

Single-use 

Endothelial marker

Dr. Terry and his team have developed a technique in which a "S" is stamped on the stromal side of the Descemet membrane^{13, 14}. This establishes a reference point to ensure the graft is correctly oriented when it is inserted into the patient's eye.



20034

Angled "S" marker

Reusable



Preparing the donor graft

Punches and trephines

Punch-assisted peripheral dissection for a clear-cut cleavage plan


As Dr. Ponzin and his team have described^{7,8}, guarded blade technology can be used to facilitate peripheral dissection and identification of the cleavage plane. Moria has developed two types of punches with a guarded blade to simplify this step: the guarded punch and the DeepWell guarded punch.

Punch-assisted final trephination for diameter

The final trephination may or may not be penetrating. In the "no-touch" technique described by Dr. Melles and his team⁶, the endothelium is drawn onto a support and then undergoes penetrating trephination. The Busin Punch, trephine blade or Hanna trephine by Moria can all be used under those conditions. Other authors have shown interest in using a guarded punch to facilitate the final peeling of the graft¹². Two guarded punches by Moria can be used for these purposes.



Reference		17212DXXX	17213DXXX	17207DXXX	17200DXXX	17150DXXX	17169	17171DXXX
Name		DeepWell Punch	DeepWell blade-holder	Guarded Punch	Busin Punch	Trephine Blade	Hanna Punch	Hanna Blade
Blade design	Double-bevel blade for a clean cut	✓	✓	✓	✓		✓	✓
	360° Blade for uniform cutting	✓	✓	✓	✓	✓	✓	✓
	Guarded blade with length adapted to the depth of the well allowing cutting under the Descemet to be performed	✓	✓	✓				
	Penetrating blade				✓	✓	✓	✓
	Blade holder sold separately		✓					
	Available sizes (mm)	7.5, 7.75, 8, 8.25, 8.5, 9.5, 10	9.5, 10	7.5, 7.75, 8, 8.5, 9.5, 10	6, 6.5, 6.75, 7, 7.25, 7.5, 7.75, 8, 8.25, 8.5, 8.75, 9, 9.5, 10	6.5, 7, 7.25, 7.5, 7.75, 8, 8.25, 8.5, 8.75, 9		7, 7.25, 7.50, 7.75, 8, 8.25, 8.50, 8.75, 9, 9.5, 10, 10.5
Well design	Deep and enveloping well, covering cornea	✓						
	Lowered base for optimized working comfort	✓						
	Wide and stable base with imprint	✓						
	Suction system	Double suction		21 suction holes	21 suction holes			
	4 non-aspirating holes to facilitate the « S » marking	✓		✓	✓			
	Graft centering system	8.5mm centering groove		4 cardinal holes	4 cardinal holes			




Preparing the donor graft

2) Hydrodissection and pneumatic dissection

Pneumatic dissection, also called "reverse big bubble technique" by Dr. Zarei Ghanavati et al.¹⁵, and hydrodissection are among the core techniques developed to prepare DMEK grafts. Principle of these techniques involves using a cannula to lift off the Descemet membrane by air (pneumatic dissection) or liquid (hydrodissection). The cannula can be inserted at different locations depending on whether the operator is performing hydrodissection^{16, 17} or pneumatic dissection^{15, 18}.

As with manual dissection, there are variations on those two techniques. For example, Dr. Ponzin et al.¹⁶ has developed a hydrodissection technique called "SubHys", which requires use of an artificial chamber and a trephine.



7504

Rycroft cannula

35G orifice, external diameter of 30G
For injection of air or liquid
Reusable



18153

Hydrodissection cannula

31G orifice, external diameter of 25G, Flat, blunt tip
Reusable



19161

Artificial chamber - Base

To maintain the donor's cornea
Reusable



19162

Artificial chamber - Cover

Compatible with base of artificial chamber (19161)
Reusable



19182

Artificial chamber

To maintain the donor's cornea
Single-use ⓘ



17204

Artificial chamber for single-use trephines

To maintain the donor's cornea
Compatible with 17201DXX and 17202DXX trephines
Single-use ⓘ



17201DXXX

Simple trephine

Available from 7mm to 9mm (increment 0.25mm) and 9.5mm
Compatible with artificial chamber (17204)
Single-use ⓘ



17202DXXX

Adjustable trephine

Pre-setting of the desired depth
Available in 6mm, 6.5mm to 9mm (increment 0.25mm), 9.5mm, 10mm
Compatible with artificial chamber (17204)
Single-use ⓘ



17150DXXX


Trephine blade

For penetrating trepanation
Available in 6.5mm, 7mm - 9mm (increment 0.25mm)
Reusable

Preparing the donor graft

3) Muraine technique

Prof. Muraine's technique consists in trephining the endothelium incompletely with two opposite hinges using a dedicated device: the Muraine punch by Moria. Donor tissue is then mounted on an artificial chamber, the endothelium upward. A Rycroft cannula is then used to hydrodissect the Descemet membrane¹⁸. This technique is described step by step in our brochure #65057 available on our website (www.moria-surgical.com).



Shape of the incision



17209D800

Muraine Punch

Allows a partial trephination with two opposite hinges

Single-use ⓘ



19161

Artificial chamber - Base

To maintain the donor's cornea

Reusable



19162

Artificial chamber - Cover

Compatible with base of artificial chamber (19161)

Reusable




19182

Artificial chamber

To maintain the donor's cornea

Single-use ⓘ



7835

Bonn forceps

0.1mm Micro-teeth

4mm Platforms

Reusable




9605

Troutman curved forceps

Active part with 7.5 mm platform.

Reusable





17225X10

Troutman forceps

Curved, 3.5mm platforms

Single-use ⓘ



7504

Rycroft cannula

35G orifice, external diameter of 30G

For injection of air or liquid
Reusable



Preparing the recipient eye

1) Keeping the patient's eye open



19078

Colibri Speculum
Lid-blades, 16mm
Reusable



20035

Adjustable speculum
Rounded lid-blades, 14mm
Reusable



18195

Schapira Speculum
15mm open lid blades
Reusable



17508x5
17508x10

Adjustable speculum
15mm lid blades
Sold in box of 5 or 10
Single-use ☒

2) Measuring and marking diameter



12994

Sourdille caliper
16mm opening
Graduation every 1mm
Reusable



19095/800
19095/850
19095

Corneal markers
Reusable

> Available in:
Diameter: 8 mm
Diameter: 8.5 mm
Diameter: 9 mm

Corneal markers
Available in boxes of 10 units
Single-use ☒

> Available in:



17518x10

Diameter: 8 mm



17519x10

Diameter: 8.5 mm



17520x10

Diameter: 9 mm



Preparing the patient's eye

3) Maintaining the patient's anterior chamber



19092

Chamber Maintainer 20G

Length, 170mm
External diameter, 1.65mm
5mm active part with internal diameter 0.90mm
Reusable

4) Descemetorhexis and removal of endothelium



19097

Gorovoy forceps

Descemetorhexis forceps
Length 11.7cm
Blunt tip
Reusable



19091

Price hook

Descemetorhexis hook
Inverted Sinskey shape
Reusable



19077/A

Spatulas

For peeling the Descemet membrane



19077/B

A: 45° angled tip
B: 90° angled tip
Reusable



19083/A

Irrigating spatulas

For peeling the Descemet membrane



19083/B

A: 45° angled tip
B: 90° angled tip
Single-use ☒



17302x5

Price hook


Descemetorhexis hook
Inverted Sinskey shape
Sold in box of 5 units
Single-use ☒



17303x5
17303x10

90° Spatula

For peeling the Descemet membrane
90°angled tip
Sold in box of 5 or 10 units
Single-use ☒



Bibliography

- 1 **Melles**, G. R., San Ong, T., Ververs, B., & van der Wees, J. (2006). Descemet membrane endothelial keratoplasty (DMEK). *Cornea*, 25(8), 987-990.
- 2 **Kruse**, F. E., Laaser, K., Cursiefen, C., Heindl, L. M., Schlötzer-Schrehardt, U., Riss, S., & Bachmann, B. O. (2011). A stepwise approach to donor preparation and insertion increases safety and outcome of Descemet membrane endothelial keratoplasty. *Cornea*, 30(5), 580-587.
- 3 **Price Jr**, F. W., & Price, M. O. (2013). Evolution of endothelial keratoplasty. *Cornea*, 32, S28-S32.
- 4 Anshu, A., Price, M. O., & **Price Jr**, F. W. (2012). Risk of corneal transplant rejection significantly reduced with Descemet's membrane endothelial keratoplasty. *Ophthalmology*, 119(3), 536-540.
- 5 Giebel AW, **Price FW** . Descemet membrane endothelial keratoplasty: the bare minimum. In: Price FW, Price MO editors. *DSEK: What You Need to Know About Endothelial Keratoplasty*. Slack Incorporated: Thorofare, NJ, 2009 pp 119–146.
- 6 Groeneveld-van Beek, E. A., Lie, J. T., van der Wees, J., Bruinsma, M., & **Melles**, G. R. (2013). Standardized 'no-touch' donor tissue preparation for DALK and DMEK: harvesting undamaged anterior and posterior transplants from the same donor cornea. *Acta Ophthalmologica*, 91(2), 145-150.
- 7 Ruzza, A., Parekh, M., Salvalaio, G., Ferrari, S., & **Ponzin**, D. (2016). Development of a new superficial punch for Descemet's Membrane Endothelial Keratoplasty donor tissue preparation. *British Journal of Ophthalmology*, 100(4), 443-445.
- 8 Parekh, M., Baruzzo, M., Favaro, E., Borroni, D., Ferrari, S., **Ponzin**, D., & Ruzza, A. (2017). Standardizing Descemet membrane endothelial keratoplasty graft preparation method in the eye bank—experience of 527 Descemet membrane endothelial keratoplasty tissues. *Cornea*, 36(12), 1458-1466.
- 9 Maharana, P. K., Sahay, P., Singhal, D., Sharma, N., & **Titiyal**, J. S. Donor preparation in descemet membrane endothelial keratoplasty. *New Front Ophthalmol*, 2019. Vol 5: 1-6.
- 10 Maharana, P. K., Sahay, P., **Titiyal**, J. S., & Sharma, N. (2019). Sinsky hook assisted roll preparation (SHARP): a modified technique for Descemet membrane endothelial keratoplasty donor preparation. *Saudi Journal of Ophthalmology*, 33(1), 28-33.
- 11 Price, M. O., & **Price**, F. W. (2013). Descemet's membrane endothelial keratoplasty surgery: update on the evidence and hurdles to acceptance. *Current Opinion in Ophthalmology*, 24(4), 329-335.
- 12 Panthier, C., Courtin, R., Gatinel, D., & Saad, A. Description de notre technique de préparation des greffons endothéliaux dans la DMEK. *Images en Ophtalmologie*, 2015. Vol. IX -n°6, 210-213
- 13 Veldman, P. B., Dye, P. K., Holiman, J. D., Mayko, Z. M., Sáles, C. S., Straiko, M. D. & **Terry**, M. A. (2015). Stamping an S on DMEK donor tissue to prevent upside-down grafts: laboratory validation and detailed preparation technique description. *Cornea*, 34(9), 1175-1178.
- 14 Veldman, P. B., Dye, P. K., Holiman, J. D., Mayko, Z. M., Sáles, C. S., Straiko, M. D. & **Terry**, M. A. (2016). The S-stamp in Descemet membrane endothelial keratoplasty safely eliminates upside-down graft implantation. *Ophthalmology*, 123(1), 161-164.
- 15 Zarei-Ghanavati, S., Khakshoor, H., & **Zarei-Ghanavati**, M. (2010). Reverse big bubble: a new technique for preparing donor tissue of Descemet membrane endothelial keratoplasty. *British Journal of Ophthalmology*, 94(8), 1110-1111.
- 16 Salvalaio, G., Parekh, M., Ruzza, A., Ferrari, S., Camposampiero, D., & **Ponzin**, D. (2014). DMEK lenticule preparation from donor corneas using a novel 'SubHyS' technique followed by anterior corneal dissection. *British Journal of Ophthalmology*, 98(8), 1120-1125.
- 17 **Szurman**, P., Januschowski, K., Rickmann, A., Damm, L. J., Boden, K. T., & Opitz, N. (2016). Novel liquid bubble dissection technique for DMEK lenticule preparation. *Graefe's Archive for Clinical and Experimental Ophthalmology*, 254(9), 1819-1823.
- 18 **Feizi**, S., & Javadi, M. A. (2015). DMEK lenticule preparation using an air dissection technique: central versus peripheral injection. *European Journal of Ophthalmology*, 26(1), 6-11.
- 19 **Muraine**, M., Gueudry, J., He, Z., Piselli, S., Lefevre, S., & Toubreau, D. (2013). Novel technique for the preparation of corneal grafts for descemet membrane endothelial keratoplasty. *American Journal of Ophthalmology*, 156(5), 851-859.

References	Manufacturer	Classes	Notified body
7835, 7850A, 13161, 13160, M1809, 20022, 20021, 6062A, 20038, 20039, 20034, 17150DXXX, 17169, 7504, 18153, 19161, 19162, 9605, 19078, 20035, 18195, 19095, 19095/800, 19095/850, 19092, 19097, 19091, 19077/A, 19077/B, 19083/A, 19083/B	Moria S.A.	I	CE : self-declaration
17504x5, 17504x10, 17221, 17521x5, 17521x10, 17212DXXX, 17213DXXX, 17207DXXX, 17200DXXX, 19182, 17204, 17209D800, 17225x10, 17508x5, 17508x10, 17518x10, 17518x5, 17519x10, 17519x5, 17520x10, 17520x5	Moria S.A.	Is	CE 0459
17171DXXX, 17201DXXX, 17202DXXX, 17302x5, 17303x5, 17303x10	Moria S.A.	Ila	CE 0459
12994	Moria S.A.	Im	CE 0459



To obtain more information

Manufacturer

MORIA SA
15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 00
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: + 1 (215) 230 7670
www.moria-surgical.com

Moria Japan K.K.
Arcadia Building 6F
1-12-3 Kanda
SudachoChiyoda-Ku
Tokyo 101-0041
JAPAN
Phone: 81-3-6260-8309
Fax: 81-3-6260-8310
www.moriajapan.com

**Moria Commercial
(CHINA) CO., LTD.**
上海目利亚贸易有限公司
Room H 6F Kaili Building NO.432
West Huai Hai Road, Changning
district, Shanghai, 200052, P.R.C
中国上海市长宁区淮海西路432号
凯利大厦6楼08室
Tel/Fax: +86 021 52586095
www.moria-surgical.com



Moria
Ophthalmic Instruments