
EN - ALTK Artificial Chamber for ONE & CBM microkeratomes USER MANUAL.....1-10	NL - Kunstmatige oogkamer ALTK voor microkeratomen ONE & CBM GEBRUIKERSHANDLEIDING51-60
FR - Chambre Artificielle ALTK pour microkératomes : ONE & CBM-ALK MANUEL UTILISATEUR..... 11-20	PT - Câmara Artificial ALTK para micro-ceratótomos ONE & CBM MANUAL DO UTILIZADOR.....61-70
ES - Cámara artificial ALTK para microqueratomes ONE & CBM MANUAL DE USUARIO21-30	SV - Artificiell Kammare ALTK för mikrokeratomerna ONE & CBM ANVÄNDARMANUAL.....71-80
IT - Camera artificiale ALTK per microcheratomi ONE & CBM MANUALE D'USO.....31-40	TR - ONE & CBM Mikrokeratomlar için ALTK Yapay Kamara KULLANICI KILAVUZU.....81-90
DE - Künstliche Kammer ALTK für Mikrokeratome ONE & CBM BENUTZERHANDBUCH.....41-50	












Ref.: 65009

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




MORIA S.A. – 15, rue GEORGES BESSE - 92160 ANTONY – France
#65009-G-01.2017

CONTENTS

I. GENERAL RECOMMENDATIONS	2
II. DESCRIPTION OF EQUIPMENT	
1. Contents (excluding console, microkeratomes & accessories)	2
2. Base of the ALTK artificial chamber (#19161)	2
3. Chamber covers for the artificial chamber	3
4. Guide rings	3
5. «ALTK» cutting lenses	3
III. TECHNICAL SPECIFICATIONS	4
IV. CHAMBER ASSEMBLY INSTRUCTIONS	
1. Preparing the «donor» graft	5
2. Adjusting the diameter of the graft cut	6
3. «Donor» graft cut with the microkeratome	6
V. TROUBLESHOOTING	7
VI. CLEANING AND STERILIZATION	
1. Disassembly	7
2. Disinfection	7
3. Double cleaning	8
4. Sterilization	8
5. Inspection and storage	9
VII. WARRANTY	
1. Scope of guarantee	9
2. Exclusions	10
3. Warranty period	10
4. Liability	10

I. GENERAL RECOMMENDATIONS

The ALTK artificial chamber has been designed for optimal use, provided the cleaning recommendations are fully respected. If, for any reason, this equipment does not work correctly, have it checked immediately by Moria. **We strongly recommend that this equipment be checked each year by Moria or its local representative.**

Any usage in the device of single-use materials and/or components other than those supplied by Moria will immediately nullify the Moria guarantee.

Before any use of the ALTK artificial chamber, you are strongly advised to familiarize yourself with the use of ONE and CBM microkeratomes, as well as the EVOLUTION console. **Please refer to each corresponding user manual (ONE microkeratome = manual #65007 for the English version and manual #65006 for the French version / CBM microkeratome = manual #65005 / Evolution console = manual #65060/INTL).**

The latest version of the user manuals and additional information about our systems are available on the Moria website <http://www.moria-surgical.com>.

II. DESCRIPTION OF EQUIPMENT

1. CONTENTS (EXCLUDING CONSOLE, MICROKERATOMES & ACCESSORIES)

Description	MORIA Reference
Base of the ALTK artificial chamber	19161
Universal chamber cover for ALTK artificial chamber («ONE» or «CBM-ALTK» type)	19162
Chamber cover for artificial chamber ALTK+ «CBM-ALTK»	19172
Guide ring for ONE microkeratome	19164
Guide ring for CBM microkeratome	19171
Tonometer	19042
ALTK cutting lens (7, 9, 11 mm)	19165
ALTK cutting lens (6, 8, 10 mm)	19166
ALTK sterilization box	22519167
ALTK cutting lens storage box	22519169

2. BASE OF THE ALTK ARTIFICIAL CHAMBER (#19161)

It consists of:

- a central piston (A)
- a graft clamping ring (B)
- a three-pin bayonet socket (C)
- a Luer connector for the infusion line (D)

Onto this common base (E) three types of artificial chamber can be mounted:

- the “ONE” type chamber, to be used with the ONE microkeratome: Assemble the universal chamber cover #19162 + «ONE» guide ring #19164.
- the “CBM-ALTK” type chamber, to be used with the CBM microkeratome:
 - Assemble either the universal chamber cover #19162 + “CBM-ALTK» guide ring #19171.
 - or the CBM-ALTK chamber cover #19172 only, as the guide is incorporated in the chamber cover.

3. CHAMBER COVERS FOR THE ARTIFICIAL CHAMBER

There are two types of design:

- **Chamber cover #19162 (F):** This can be used to make a graft cut with “ONE” and “CBM-ALTK” microkeratomes. You must use:
 - the guide ring #19164 with the “ONE” system (L),
 - the guide ring #19171 with the “CBM-ALTK” system (M).In addition, this chamber cover allows you to set the diameter of the graft cut via an adjustment ring.
- **Chamber cover #19172:** This incorporates the CBM microkeratome “slide” (G) and pivot (H) system. No guide ring is necessary. This chamber cover can only be used with the CBM microkeratome. The graft cut is automatically as large as possible. It is therefore necessary to use a trephine, as a second step, in order to adjust the diameter of the graft cut on the eye of the recipient.

Mounting chamber covers on the base of the artificial chamber (the photos provided as an example show chamber cover #19172).

The chamber covers (“ONE” and “CBM-ALTK” type) are placed on the base of the artificial chamber to hold the graft (I). Turn clockwise (J). The three-pin bayonet system allows you to lock the chamber cover in position (K).

4. GUIDE RINGS (NOT TO BE USED WITH UNIVERSAL CHAMBER COVER #19162)

There are two types of guide ring:

- # 19164: type “ONE” for use with the ONE microkeratome (L)
- # 19171: type “CBM-ALTK” for use with the CBM microkeratome (M).

The locking screw on the guide ring allows you to secure the chamber cover in place (N).

5. “ALTK” CUTTING LENSES

Two cutting lenses are available (O). They are used to check the diameter of the punch that will be made on the graft.

- Lens #19165 indicates 7mm, 9mm and 11mm diameters.
- Lens #19166 indicates 6mm, 8mm and 10mm diameters.

III. TECHNICAL SPECIFICATIONS

	Dimensions	Weight	°C/storage	RH / storage
Base of the ALTK artificial chamber #19161 + chamber cover #19172 (device assembled)	Diam = 80 mm H = 60 mm	334 g	Between 15°C/59°F and 35°C/95°F	Between 30 and 75 %

IV. CHAMBER ASSEMBLY INSTRUCTIONS

Before any usage, confirm that the device is suitable and operable for the proposed medical procedure. The device must not be subjected to any shock and must be free of marks or traces of wear, in particular the parts in contact with the patient or medical staff.

Check that the device will not interfere with other medical devices used during the procedure.

1. PREPARING THE «DONOR» GRAFT

Only the preparation of the donor graft is detailed in this manual. For the preparation of the patient eye, please refer to the manual of your equipment (#65005 for the CBM microkeratome ; #65007 (English version) or #65006 (French version) for the ONE microkeratome).

ACTIONS	PRECAUTIONS
1. Place the base of the artificial chamber on the table and make sure it is stable.	Make sure that the central piston is in the lowest position by turning the graft clamping ring counterclockwise.
2. Connect the infusion line to the Luer connector (P).	Make sure that the Luer connector is secure and airtight. The ALTK chambers are supplied with silicone tubing. We recommend replacing the silicone tubing after each surgery (reference: 117.18 for a meter of silicone tube).
3. Connect the other end of the infusion line to the infusion bottle.	
4. Fill the chamber before irrigation. Temporarily halt the infusion until the graft has been put in place and the chamber closed.	
5. Position the graft on the top of the piston, endothelium downwards.	
6. Start the irrigation and, with the irrigation flowing, move the graft to remove any air bubbles present under the cornea. Once all bubbles have been removed, stop irrigation. Center the cornea on the piston.	You must use corneas with a minimum diameter of 16 to 17mm in order to have enough sclera to ensure airtightness when the chamber is closed.
7. Position the chamber cover on the base (Q). Insert the pins of the chamber cover into the slots in the base. Then turn the chamber cover about 15° counterclockwise.	
8. Check that the graft is properly centered. If the cornea is off-center, gently remove the chamber cover, recenter the cornea and repeat the above procedure.	
9. Pull the central piston and cornea upward by turning the graft clamping ring clockwise. The central piston must be raised until the donor cornea comes into contact with the top orifice of the chamber cover (R).	The system must be firmly adjusted, solely by hand, in order to ensure that it is perfectly airtight and to hold the graft during the punch.
10. Start the irrigation.	Check that the system is airtight. If it is not, lock the graft clamping ring clockwise.
11. Arrange the infusion bottle so as to obtain a pressure of at least 65 mmHg in the artificial chamber . Check the pressure using a tonometer (S). Adjust the height of the bottle if necessary.	The cornea surface in contact with the tonometer must be contained in the area determined by the inscribed circle.

2. ADJUSTING THE DIAMETER OF THE GRAFT CUT (ON CHAMBERS “ALTK+ONE” AND “ALTK+CBM-ALTK”)

ACTIONS	PRECAUTIONS
1. Position the adjustment ring on the chamber cover, with the marking arrows “+” and “-” (T) upwards, and screw it down fully.	Be careful to tighten only finger-tight.
2. Position the desired cutting lens (U).	The lenses supplied with the ONE microkeratome can also be used instead of ALTK lenses. To do so, you must put in place the ONE guide ring #19164. The lenses supplied with the CBM microkeratome can also be used instead of ALTK lenses. To do so, you must put in place the CBM-ALTK guide ring #19171.
3. Check the cut surface. The cut region corresponds to the internal limit of the reticle.	
4. Adjust the cut surface to the desired diameter. By turning the adjustment ring clockwise, the cut surface increases (“+” arrow on the top face of the ring). By turning the adjustment ring counterclockwise, the cut surface decreases (“-” arrow on the top face of the ring). When you have obtained the desired diameter, remove the cutting lens.	 WARNING DO NOT CHANGE THE POSITION OF THE ADJUSTMENT RING AFTERWARD!
5. Position the guide ring on the chamber cover (V example of an «ALTK+ONE» chamber cover). Keep it in position by pressing down on it and secure it using the locking screw. The system is now ready for cutting with the microkeratome.	Make sure that the guide ring and the chamber cover are properly <u>in contact and parallel</u> in order to obtain the correct cut thickness and diameter.  If the guide ring is incorrectly positioned, the thickness and diameter of the cut can change.

3. “DONOR” GRAFT CUT WITH THE MICROKERATOME

ACTIONS	PRECAUTIONS
1. Insert the microkeratome into the slide and begin to cut. Advance it using a steady movement and maintain the movement until the microkeratome comes exits the slide on the other side (W: chamber ALTK+ONE).	<u>For chamber ALTK+ONE:</u> make sure that the stop is removed from the «ONE» microkeratome, as you want to obtain a complete resection. <u>For chamber ALTK+CBM-ALTK:</u> The heads of the CBM-ALTK microkeratome do not have any hinge stops so as to obtain a complete resection. These heads must under no circumstances be used for a LASIK procedure.
2. Remove the cutting lens from the head of the microkeratome using small toothless pliers (X: chamber ALTK+ONE).	

Note: Chamber cover version CBM-ALTK #19172 (without guide ring)

This version allows you to cut a maximum diameter graft without first needing to adjust it. After positioning the chamber cover on the base, the procedure can begin directly with the graft cut. Insert the CBM-ALTK head into the slide and begin to cut. Advance it using a steady movement and maintain the movement until the CBM-ALTK head exits the slide on the other side. You must obtain a complete resection. Gently remove the cutting lens from the CBM-ALTK head using small toothless pliers.

V. TROUBLESHOOTING

- If the device does not function properly, check for any obvious faults.
- If there is a leaking during irrigation, check the state of the tubing and replace it if necessary.
- Check that the tubing is properly in place, from the Luer connector to the end of the tube inserted in the elbow under the chamber.
- If there is a problem securing the guide ring on the chamber cover, check the state of the screw threads.
- If the guide rings stick in the slides, check for any deposit or obstruction in these slides.

VI. CLEANING AND STERILIZATION

Recommendations for cleaning, disinfection and sterilization are updated and available on the Moria website: <http://www.moria-surgical.com>.

For cleaning and sterilization recommendations for the console, microkeratomes, and their accessories, please see the corresponding user guides.

WARNING: The following information is provided as recommendations. They must be adapted to suit the laws in force in the country of use. Products, in particular for cleaning and disinfecting, must be used according to the product manufacturer's instructions.

1. DISASSEMBLY

All parts of the device must be disassembled before disinfecting and cleaning (the chamber, chamber cover and guide ring must be treated separately).

Please refer to the instructions below for how to **disassemble the artificial chamber**:

- Unscrew the guide ring locking screw (1/2 turn) (Y) and pull out the guide ring (Z).
- Unscrew the diameter adjustment ring from the chamber (AA & AB).
- Turn the graft clamping ring counterclockwise until you feel resistance (AC). Do not force it.
- While you are doing this, you can see at the top of the chamber the central piston descending. When you feel resistance, the central piston has reached its bottom position. Move the piston back to middle position so as not to block the system.
- Turn the chamber cover about 15° until you can free the bayonet and lift off the chamber cover (AD & AE). The chamber is now disassembled. The other parts that make up the base should not be disassembled.

2. DISINFECTION

All parts of the device must be immersed immediately upon being taken out of the block in a non-aldehyde detergent-disinfectant solution for a minimum of 15 minutes, in accordance with the manufacturer's instructions.

The base of the artificial chamber must be ventilated and rinsed several times with decontaminant fluid and make sure that no foreign substances remain.

3. DOUBLE CLEANING

- 1st manual cleaning * (or use an automatic washer, provided the machine does not recycle the cleaning fluid).
- Carefully rinse using tap water.
- Immerse the device in a non-aldehyde detergent-disinfectant solution (different or prior bath).
- 2nd manual cleaning * (or use an automatic washer, provided the machine does not recycle the cleaning fluid).
- Carefully rinse with filtered distilled or osmosis water.
- Wipe the device with a non-fabric, lint-free, single-use swab. Once the components are completely dry, package them promptly to avoid contamination. Dry the inside of the base of the device using filtered compressed air.

Important: Change the bath after each use.

** Manual cleaning: Mechanically clean the device for a minimum 2 minutes using soft brushes. (You must not use metal brushes). Do not use abrasives or scrapers to clean the parts. Doing so will reduce precision, produce scratches, and damage the metal surfaces.*

Note about unconventional transmissible agents (UCTAs):

For France, follow the recommendation of the Health Ministry concerning the deactivation of UCTAs (DGS/5 C/DHOS/E 2 n° 2001-138 of 14/03/01).

- Immerse the device immediately after removing from the block into a non-aldehyde preferably enzymatic detergent-disinfectant solution, for a minimum of 15 minutes, according to the manufacturer's recommendations, then rise the devices.
- Then immerse the device in a soda solution (1 N) for one hour. **Only the tonometer must not be cleaned with soda.**
- Finally, rinse manually three times in succession (with tap water), with a pH check at the first rinse to avoid any risk of burning (pH neutral). The final rinse must be done with microbiologically controlled water.
If the pH is not neutral, rinse several more times until the pH is acceptable (close to pH7).
- Wipe the device with a non-fabric, lint-free, single-use swab and put it back into the normal disinfection/sterilization circuit.

4. STERILIZATION

Sterilize (**except for cutting lenses**) using a self-regulating sterilizer (steam sterilizer) at minimum 134°C (273°F) for 18 minutes.

To remove the water, it is very important to use the self-regulating sterilizer's drying program.

The parts must be removed from the self-regulating sterilizer immediately after sterilization and must be cooled to room temperature before use.

Avoid thermal shock, as it can weaken the material.

IMPORTANT ! Leaving residual humidity in the ALTK chamber can produce oxidation. Make sure that the parts are completely dry before storing them, by passing them through the self-regulating sterilizer's drying program (sterilizers without a drying cycle are not recommended).

5. INSPECTION AND STORAGE

It is important to check that the devices are completely clean, dry, and free of any residue (organic or otherwise). Stained devices must immediately be removed from instrument boxes and must not be used in that state.

All parts must then be stored in their dedicated boxes, in an appropriate, dry place.

Do not store the instruments near products or an environment that may be corrosive or magnetic. Prevent all contact between the instruments particularly if they are made of different materials.

VII. WARRANTY

1. SCOPE OF WARRANTY

For the console, the microkeratomes, and their accessories, please refer to the corresponding user manuals.

Description	MORIA Reference
Base of the ALTK artificial chamber	19161
Universal chamber cover for ALTK artificial chamber («ONE» or «CBM-ALTK» type)	19162
Chamber cover for artificial chamber ALTK+ «CBM-ALTK»	19172
Guide ring for ONE microkeratome	19164
Guide ring for CBM microkeratome	19171
Tonometer	19042
ALTK cutter lens (7, 9, 11mm)	19165
ALTK cutter lens (6, 8, 10mm)	19166
ALTK sterilization box	22519167
ALTK cutter lens storage box	22519169

- The warranty covers spare parts and labor necessary to bring the equipment to proper operating condition. Any equipment to be returned **must first be disinfected**, then sent in its original packing.
- Maintenance and replacement of spare parts must be done exclusively by MORIA-approved staff.
- Any defective part replaced during the warranty period remains the property of MORIA.
- The warranty excludes non-reusable instruments and accessories (such as single-use tubing).

2. EXCLUSIONS

- Defect or malfunction in the system appearing outside the warranty period (defined in the “warranty period” paragraph).
- Normal wear and tear.
- Negligence, user fault, or any use that does not conform to user-guide specifications.
- Console connected to an inappropriate energy source.
- The use of consumables, spare parts or accessories other than those supplied by MORIA.
- Any disassembly, modification, or intervention done on the equipment by a person not authorized by MORIA.

3. WARRANTY PERIOD

- The warranty takes effect on the day the equipment is shipped.
- The warranty is for a period of 12 months from the effective date.
- In addition to the contractual one-year warranty, repairs are guaranteed for three months from the date that they are invoiced.

4. LIABILITY

- MORIA's liability is limited to the supply of services set out in the «scope of warranty» paragraph. MORIA shall not be held liable for any direct or indirect damage to a client caused by interventions under this warranty.
- Any dispute over the interpretation or execution of this contract or of these general terms and conditions shall be governed by the Commercial Court of Nanterre (France).

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