CE 0459









Ref.: 65060/INTL

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










CONTENTS

I. DISCLAIMER	2
A. MAINTENANCE AND WARRANTY	2
B. USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES	2
II. GUIDANCE AND MANUFACTURER'S DECLARATION: ELECTROMAGNETIC EMISSIONS AND IMMUNITY	2
III. EQUIPMENT AND ACCESSORIES LIST	2
A. EQUIPMENT LIST	2
B. ACCESSORIES LIST	3
IV. TECHNICAL SPECIFICATIONS	6
V. REGULATORY INFORMATION	6
VI. INTENDED USE	6
VII. LABELING INFORMATION	7
VIII. INSTALLATION AND CONNECTION	7
A. INSTALLATION RECOMMENDATIONS	7
B. ELECTRICAL SETTING	7
C. BATTERY TEST	8
D. CONNECTING THE ASPIRATION TUBING	8
E. VACUUM TEST	8
F. LOW VACUUM FUNCTION	9
G. MANUAL MODE (CB, CBm, CBSU, ONE, ONE USE)	9
1. PNEUMATIC TURBINE CONNECTIONS	9
2. DISASSEMBLY OF THE NITROGEN/AIR SUPPLY	10
H. AUTOMATED MODE (M2, M2SU, ONE USE-PLUS, EPI-K™)	10
IX. USE OF FOOTPEDAL	11
A. COMPATIBILITY FOOTPEDAL-CONSOLE	11
B. COMPATIBILITY FOOTPEDAL-HANDPIECE	11
C. MANUAL MODE	11
D. AUTOMATED MODE	11
X. TROUBLE SHOOTING	11
A. POWER SUPPLY	11
B. VACUUM LOSS	12
C. RESIDUAL VACUUM	12
D. TURBINE PRESSURE (MANUAL MODE)	12
XI. CARE & MAINTENANCE	12
A. HOW TO SELECT THE OPERATING VOLTAGE (115V OR 230V POSITION)	12
B. INSTALLING A NEW BATTERY	12
C. REPLACING FUSES	13
XII. WARRANTY	13
A. SCOPE OF WARRANTY	13
B. NON-APPLICATION OF WARRANTY	13
C. WARRANTY PERIOD	13
D. LIABILITY	13

XIII. DRAWINGS	13
A. CONTROL UNIT (FRONT)	13
B. CONTROL UNIT (BACK)	13
C. FOOTPEDAL	13
D. CONNECTING AND DISCONNECTING NITROGEN/AIR HOSE (#19120)	13
E. HOW TO SELECT THE OPERATING VOLTAGE (115V OR 230V POSITION)	13
F. CHANGING THE BATTERY	13
G. CHANGING FUSES	13
H. INSTALLING TUBING	13

This instruction manual is for EVOLUTION 3E with serial numbers 5000 and above.
For the EVOLUTION 3E with serial numbers below this number, please refer to user manual (#65051).
The most recent version of this user guide and additional information on your keratome are available on MORIA website: <http://www.moria-surgical.com>.

I. DISCLAIMER

A. MAINTENANCE AND WARRANTY

The EVOLUTION 3E system has been designed for optimal operation, provided that the recommendations listed in this user manual are followed carefully.
The lifetime of the Evolution 3E system is 5 years. This lifetime can be extended to 15 years if the maintenance plan described above is strictly followed. The maintenance operations in the fifth and tenth years (M5 and M10) must be carried out by Moria. They guarantee the electrical safety and operating performance of the Evolution 3E Console and thus the user's and the patient's safety.
The other annual, intermediate maintenance operations will be carried out by Moria or by its representative who will have been previously trained and authorized by Moria. This representative must have a valid training certificate and must use only spare parts supplied by Moria.
MORIA shall not guarantee the performance of the Evolution 3e system and patient and user safety if the maintenance plan is not followed.

Years	Ref.	Check	Electrical safety	Replacement of the primary circuit	Battery change	Turbine gasket change	Pump change	Tubing change	Engine board change	Vacuum solenoid valve change
1	M1	yes	no	no	no	yes	no	yes	no	no
2	M2	yes	no	no	no	yes	no	no	no	no
3	M3	yes	no	no	no	yes	no	yes	no	no
4	M4	yes	no	no	no	yes	no	no	no	no
5	M5	yes	yes	no	yes	yes	yes	yes	no	no
6	M6	yes	no	no	no	yes	no	no	no	no
7	M7	yes	no	no	no	yes	no	yes	no	no
8	M8	yes	no	no	no	yes	no	no	no	no
9	M9	yes	no	no	no	yes	no	yes	no	no
10	M10	yes	yes	yes	yes	yes	yes	yes	yes	yes
11	M11	yes	no	no	no	yes	no	no	no	no
12	M12	yes	no	no	no	yes	no	yes	no	no
13	M13	yes	no	no	no	yes	no	no	no	no
14	M14	yes	no	no	no	yes	no	yes	no	no
15	No maintenance									

If, for any reasons, the system does not perform properly, have it checked immediately by MORIA. MORIA strongly recommends having the system thoroughly inspected by MORIA on a routine basis every year.

The use of materials and/or components of a brand other than MORIA with the EVOLUTION 3E system will immediately nullify the MORIA warranty. MORIA may not be held responsible for any damage resulting from the use of materials and/or components of a brand other than MORIA.

As only MORIA and its agents are fully expert in MORIA products, servicing and maintenance must be carried out by MORIA or its approved agents.

MORIA shall not be held liable for any malfunction or damage to the apparatus, poor results, or surgical complications due to maintenance being having been carried out by an unqualified operator or third party.

Any such unauthorised intervention shall render the guarantee and any maintenance contract null and void.

B. USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES

Single-use devices should not be re-used. Doing so will negatively affect their clinical performance and increase the potential for adverse events.
The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and dama ge the microkeratome.

MORIA shall not be held liable in the event of a malfunction or damage to the microkeratome, poor results or surgical complications due to the reuse of a single-use product, or the use of consumables other than those supplied by MORIA.

MORIA handpieces must only be connected to MORIA devices (console unit, heads, suction rings, etc.).
All warranties become null and void if the microkeratome degrades or malfunctions due to such practices.
The most recent version of this user guide and additional information on your keratome are available on MORIA website: <http://www.moria-surgical.com>.

II. GUIDANCE AND MANUFACTURER’S DECLARATION: ELECTROMAGNETIC EMISSIONS AND IMMUNITY

Refer to annexe document (#65073).

III. EQUIPMENT AND ACCESSORIES LIST

A. EQUIPMENT LIST

Designation	Regulatory designation	MORIA reference
EVOLUTION 3E Console (S/N above 5000)	EVOLUTION 3e CONTROL UNIT	19380
EVOLUTION 3-3E Control Footpedal	FOOT PEDAL FOR EVOLUTION	19361
EVOLUTION 3E Control Footpedal Epi-K™	EVOLUTION 3e FOOTSWITCH	19381
Evolution 3E Footswitch - China	EVOLUTION 3e FOOTSWITCH	19381C
Evolution 3E footswitch - Japan	EVOLUTION 3e FOOTSWITCH	19381J
EVOLUTION 3 Supply Cords (CEE) (2.50m) / Cable (CEE)	EVOLUTION 3 SUPPLY CORD (EEC)	19362
EVOLUTION 3 Supply Cords (USA) (2.50m) / Cable (USA)	EVOLUTION SUPPLY CORD (USA)	19363
EVOLUTION 3 Supply Cords (UK) (2.50m) / Cable (UK)	EVOLUTION 3 SUPPLY CORD (UK)	19364
EVOLUTION 3 Supply Cords (China) (2.50m) / Cable (China)	EVOLUTION 3 SUPPLY CORD CHINA	19516
EVOLUTION 3 Supply Cords (Brazil) (2.50m) / Cable (Brazil)	EVOLUTION 3 SUPPLY CORD BRAZIL	19521
Supply cord USA	SUPPLY CORD (USA)	19451
Carrying Case	N/A	19511
User manual	N/A	65060/INTL
User manual - China	N/A	65060/ZH
User manual- Brazil	N/A	65060/BR
Annexe "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	N/A	65073

B. ACCESSORIES LIST

Designation	Regulatory designation	MORIA reference	User manual	
Dry medical grade nitrogen gaz / Dry and filtered medical grade air Gas Supply Hose* *This hose is supplied when the EVOLUTION 3E unit is ordered together with MORIA manual microkeratomes"	GAS SUPPLY TUBING	19120		
Tonometer	BARRAQUER autoclavable tonometer	19042		
Aspiration tubing	Evolution 2-3-3E tubing x10	19138		
Turbine hose (for ONE or CBm-ALTK turbine)	SHORT CONNECTOR TURB FLEX HOSE	19353		
Turbine handpiece for ONE microkeratome	GAS TURBINE FOR LSK ONE	19155		
Turbine handpiece for CB microkeratome	LSK C & B TURBINE MOTOR	19303		
M2 handpiece	M2 ELECTRIC MOTOR	19326		
Epi-K™ handpiece (green color)	EPI-K HANDPIECE	19342		
One Use- <i>Plus</i> handpiece (blue color)	ONE USE PLUS MOTOR	19345		
One Use- <i>Plus</i> One-Handed handpiece (blue color)	ONE HAND ONE USE PLUS HANDPIECE	19345OH		
One Use- <i>Plus</i> hanpiece for artificial chamber (black color)	OUP HANDPIECE F/ ARTIFCIAL CHAMBER	19175		
Reusable Epi-K™ metal ring size -1, 0, +1 with stops at 7.5 - 8.0 - 8.5 - 9.0 (green color)	EPI-K REUSABLE SUCTION RING -1	19389/-1	65043	
	EPI-K REUSABLE SUCTION RING 0	19389/0		
	EPI-K REUSABLE SUCTION RING +1	19389/1		
Reusable Epi-K™ metal ring size Large Diameter (LD) -1 with stops at 8.0 - 8.5 - 9.0 - 9.5 (green color)	EPI-K LD REUSABLE SUCTION RING -1	19399/-1		
Box of 10 sterile Epi-K™ heads	EPI-K HEAD	19390		
Box of 10 Epi-K™ ring -1 and 10 Epi-K™ heads	EPI-K SET (-1 RING)	19384	65044	
One Use- <i>Plus</i> : ring –1 and Head SBK (90) or 130	ONE USE PLUS SBK SET (90/-1)	19336/90	65039	
	ONE USE PLUS SET (130/-1)	19336/130		
One Use- <i>Plus</i> : ring 0 and Head SBK (90) or 130	ONE USE PLUS SBK SET (90/0)	19337/90		
	ONE USE PLUS SET (130/0)	19337/130		
One Use- <i>Plus</i> : ring –1 Large Cut and Head Large Cut 110L and 130L	ONE USE PLUS 110 LARGE CUT SET	19354/110		
	ONE USE PLUS 130 LARGE CUT SET	19354/130		
One Use- <i>Plus</i> Large cut reusable metal ring size -1, 0 with stops at 7.5 - 8.0 - 8.5 - 9.0 (mauve color)	OUP LC -1 REUSABLE SUCTION RING	19519/-1		
One Use- <i>Plus</i> reusable metal ring size -1, 0 with stops at 7.5 - 8.0 - 8.5 - 9.0 (blue color)	OUP REUSABLE SUCTION RING -1	19391/-1	65040	
	OUP REUSABLE SUCTION RING 0	19391/0		
One Use- <i>Plus</i> reusable metal ring size +1, +2, +3 with stops at 7.0 - 7.5 - 8.0 - 8.5 (yellow color)	OUP REUSABLE SUCTION RING 1	19391/1		
	OUP REUSABLE SUCTION RING 2	19391/2		
	OUP REUSABLE SUCTION RING 3	19391/3		
One Use- <i>Plus</i> reusable metal oval rings size +1, +2, +3, with stops at 7.0 - 7.5 – 8.0 - 8.5 (grey color)	OVAL FLAP OUP METAL SUCTION RING +1	19391/1/OV		
	OVAL FLAP OUP METAL SUCTION RING +2	19391/2/OV		
	OVAL FLAP OUP METAL SUCTION RING +3	19391/3/OV		
Box of 10 sterile heads One Use- <i>Plus</i> SBK (90) or 130	ONE USE PLUS SBK HEAD	19393/90		
	ONE USE PLUS 130 HEAD	19393/130		
M2-M2SU suction ring (size -1, 0, +1, +2, +3)	M2 SUCTION RING -1	19325/-1		65013
	M2 SUCTION RING 0	19325/0		
	M2 SUCTION RING +1	19325/1		
	M2 SUCTION RING +2	19325/2		
	M2 SUCTION RING +3	19325/3		
M2-M2SU suction ring size -1 Large-Cut	M2 SUCTION LARGE CUT RING -1	19379/-1		
M2 metallic head	M2 HEAD 110	19327/110		
	M2 HEAD 130	19327/130		
Box of 10 sterile disposable M2 blades	M2 BLADE	19329		
	M2 BLADE X100 (CHINA)	19329C4		
	M2 BLADE X50	19329L3		
Box of 10 sterile disposable M2SU heads	M2SU 90 HEAD	19334/90		
		19334/130		
	M2SU 90 HEAD X80 (CHINA)	19334/90C5		

Designation	Regulatory designation	MORIA reference	User manual
Microkeratome CB manual suction rings (size -1, 0, +1, +2)	LSK C & B MANUAL FIXATION RING -1	19309/-1	65005
	LSK C & B MANUAL FIXATION RING 0	19309/0	
	LSK C & B MANUAL FIXATION RING 1	19309/1	
	LSK C & B MANUAL FIXATION RING 2	19309/2	
Box of 10 sterile CBm blades	CBM BLADE	19333	
	CBM BLADE X50	19333L3	
Box of 10 sterile CBSU heads	CBSU 110 HEAD	19312/110	65006 - 65007
	CBSU 130 HEAD	19312/130	
Box of 10 ONE Microkeratome blades	ONE MICROKERATOME BLADE	19137	65078
CBm-ALTK Head	CBM / ALTK HEAD 110	19170/110	
	CBM / ALTK HEAD 130	19170/130	
	CBM / ALTK HEAD 150	19170/150	
	CBM / ALTK HEAD 200	19170/200	
	CBM / ALTK HEAD 250	19170/250	
	CBM / ALTK HEAD 300	19170/300	
	CBM / ALTK HEAD 350	19170/350	
	CBM / ALTK HEAD 400	19170/400	
Box of 5 sterile CBSU heads for artificial chamber	CBSU 50 HEAD - ART. CHAMBER	19178/50	
	CBSU 90 HEAD - ART. CHAMBER	19178/90	
	CBSU 110 HEAD - ART. CHAMBER	19178/110	
	CBSU 130 HEAD - ART. CHAMBER	19178/130	
	CBSU 200 HEAD - ART. CHAMBER	19178/200	
	CBSU 250 HEAD - ART. CHAMBER	19178/250	
	CBSU 300 HEAD - ART. CHAMBER	19178/300	
	CBSU 350 HEAD - ART. CHAMBER	19178/350	
Disposable CBSU guide ring for disposable artificial chamber (19182)	CBSU GUIDE RING FOR ARTIFICIAL CHAMBER	19180	N/A
Box of 5 sterile One Use Large-Cut heads for artificial chamber	OU LC 110 HEAD ARTIFICIAL CHAMBER	19184/110	65009
	OU LC 130 HEAD ARTIFICIAL CHAMBER	19184/130	
	OU LC 150 HEAD ARTIFICIAL CHAMBER	19184/150	
	OU LC 200 HEAD ARTIFICIAL CHAMBER	19184/200	
	OU LC 250 HEAD ARTIFICIAL CHAMBER	19184/250	
	OU LC 300 HEAD ARTIFICIAL CHAMBER	19184/300	
	OU LC 350 HEAD ARTIFICIAL CHAMBER	19184/350	
	OU LC 400 HEAD ARTIFICIAL CHAMBER	19184/400	
	OU LC 450 HEAD ARTIFICIAL CHAMBER	19184/450	
	OU LC 500 HEAD ARTIFICIAL CHAMBER	19184/500	
One Use Large-Cut guide ring for ALTK artificial chamber (19161-19162)	OU LC 550 HEAD ARTIFICIAL CHAMBER	19184/550	N/A
	OU LC RING FOR ALTK ARTIFICIAL CHAMBER	19168	
	OUP LC RING F/ ALTK ARTIFICIAL CHAMBER	19173	
	OU LC GUIDE RING FOR ARTIFICIAL CHAMBER	19186	
	OUP LC GUIDE RING F/ ART. CHAMBER X5	19187	
One Use- <i>Plus</i> Large-Cut guide ring for ALTK artificial chamber (19161-19162)			
Disposable One Use Large-Cut guide ring for disposable artificial chamber (19182)			
Disposable One Use- <i>Plus</i> Large-Cut guide ring for disposable artificial chamber (19182)			

Designation	Regulatory designation	MORIA reference	User manual
Base of the ALTK artificial chamber	Donor artificial chamber base	19161	65009
ALTK/CBm cover with guide ring	ALTK/CBm cover with guide ring	19172	
Universal chamber cover for ALTK artificial chamber	Artificial chamber cover for ALTK	19162	
Infusion tubing set for artificial chamber	INFUSION TUBING SET FOR ARTIFICIAL CHAMBER	19181	N/A
Disposable artificial chamber	ONE - SU ART. CHAMBER x1	19182	
OUP Ring sterilization box	OUP ring sterilization box	22519513	65040
Polymer box for OUP handpiece	Polymere box for OUP handpiece	22519514	65039-65040
Sterilization box for M2	Sterilization box for M2	22519330	65013-65062
ALTK sterilization box	ALTK sterilization box	22519167	N/A
ALTK lenses storage box	ALTK lenses storage box	22519169	
Storage box for M2	Storage box for M2	22519331	
Storage box F/art. Chamber OUP handpiece	STORAGE BOX F/ART. CHAMBER OUP HANDPIECE	22519174	
Sterilization box for LSK	Sterilization box for LSK	22519139	
LSK CB sterilization tray	LSK CB sterilization tray	22519343	
OUP motor loosening key	OUP MOTOR LOOSENING KEY	19345C	65005-65078
Applanator lense 8-10 F/LSK C&B	Applanator lense 8-10 F/LSK C&B	19311	
ALTK applanator lense 6-8-10	ALTK applanator lense 6-8-10	19166	65009-65078
ALTK applanator lense 7-9-11	ALTK applanator lense 7-9-11	19165	
Artificial Chamber cover for ALTK/CBm	Artificial Chamber cover for ALTK/CBm	19171	
Artificial chamber guide ring for ALTK	Artificial chamber guide ring for ALTK	19164	65009
ONE head blade push	ONE head blade push	19135	N/A

IV TECHNICAL SPECIFICATIONS

The control unit is mainly composed of the following elements:

- 1 x 12V-7Ah battery (back-up)
- 1 x transformer
- 2 x vacuum pumps: vacuum circuit (main pump + « back-up » pump)
- valves for the vacuum & gas circuit
- 1 x gas regulator used for the set up of the nitrogen pressure going to the turbine (manual microkeratomes)
- Electronic boards.

Reference	mm	kg	Pressure	Elect. Supply
19380	430 x 240 x 190	13.6	IN : 4-6 bars / 58-87 psi OUT : 2.5-3.5 bars / 36-50 psi Precision of displayed values: +/-0,2bar = 2.9 psi. Vacuum display: +/- 30 mmHg	100-120V / 220-240V 50/60Hz Fuse : 500 mA high switching power Battery : 12 V - 7 Ah (Pb) ¹ Power input: 50 VA
19381	Double footswitches: 220 x 200 x 60 Single footswitch : 100 x 70 x 30	2.55		
19361	Double footswitches: 160 x 130 x 70 Single footswitch : 100 x 70 x 30	1.05		
19381C	Double footswitches: 160X130X70 Simple pédale: 100X70X30	1.05		
19381J	Double pédale: 220X200X30 Single footswitch: 100X70X30	2.55		

¹ If a battery change is necessary, the used battery is classified as waste, which must be dealt with according to the specific regulations for collection, treatment, recycling and disposal set out in directive 2006/66/EC of 9 September 2006 relating to batteries and accumulators and waste batteries and accumulators. As such the battery must be disposed of at an appropriate collection point.

The Evolution 3E control unit and the electrical motors / turbines must follow the following conditions:

	Temperature	Humidity rate (without condensation)	Atmospheric pressure
Transportation	5 - 45°C / 41°F – 113°F	30 - 90%	
Storage	10 - 40°C / 50°F – 104°F	30 - 75%	
Normal conditions of use	10 - 40°C / 50°F – 104°F	30 - 75%	70 - 106 kPa


V. REGULATORY INFORMATION

	MORIA S.A. 15, rue Georges Besse - 92160 Antony - France Phone +33 (0) 146 744 674 Fax +33 (0) 146 744 670 moria@moria-int.com http://www.moria-surgical.com
Electrical Safety Standard, electromagnetic compatibility	IEC 60601-1 Classe I type BF
Annexe "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	Refer to annexe document (#65073).
Customer Service Information	Contact your local dealer or MORIA.
EUROPE	CE 0459 according to MDD 93/42/CEE
USA	Product registered at the Food and Drug Administration (FDA). Caution USA only: Federal law restricts the use of this device to physicians or licensed practitioners.
	For EU customers only: this symbol indicates that within the European Union, the product must be discarded in a separate collection bin at the end of its useful life. This applies not only to this device, but also all accessories, including the footpedal and electrical motors, regardless of whether those accessories are marked with the symbol. Do not dispose as unsorted municipal waste. For users out of EU: please refer to local environmental regulatory regarding waste of electrical and electronic equipment.

VI. INTENDED USE

The power unit has been designed to operate Moria keratomes by means of electric motor or turbine. Evolution 3E console must be used only by experienced ophthalmic surgeons with specific knowledge and training in corneal surgery.

VII. LABELING INFORMATION

 XXXXXX	CATALOGUE NUMBER OR SKU
 XXXXXXXX	SERIAL NUMBER
 YYYY/MM/DD	DATE OF MANUFACTURE
 XXXXXXX	MANUFACTURER
 XXXXXX	CAUTION : CONSULT INSTRUCTIONS FOR USE FOR ANY IMPORTANT INFORMATION ABOUT THE SECURITY
 XXXXXX	CONSULT OPERATING INSTRUCTIONS
	DISCARD IN A SEPARATE COLLECTION BIN
	WARNING AGAINST ELECTROSTATIC DISCHARGES
	WARNING AGAINST INTERFERENCE RISK
	REFER TO USER MANUAL
	SIGN OF GENERAL SAFETY
	ALTERNATIVE CURRENT
	TYPE « BF » APPLIED PARTS: PROTECTION AGAINST ELECTRIC SHOCK : • LEAKAGE CURRENT • GROUNDED POWER SYSTEM
	NORMAL CONDITIONS OF USE EVO3E CONTROL UNIT AND MOTORS/TURBINES: TEMPERATURE: 10°C – 40°C / 50°F – 104°F, HUMIDITY RATE: 30 - 75%
	NORMAL CONDITIONS OF USE SINGLE USE ITEM : TEMPERATURE: 15°C – 35°C / 59°F – 95°F, HUMIDITY RATE: 30 - 70%
IP X6	FOOTSWITCH PROTECTED AGAINST POWERFUL WATER JETS
IP X8	FOOTSWITCH PROTECTED AGAINST WATER SUBMERSION BEYOND 1 METER (3,3 FEET)
	THIS SIDE UP
	FRAGILE
	KEEP DRY
Rx only	CAUTION FOR USA ONLY: US FEDERAL LAWS RESTRICT THIS DEVICE TO SALE BY, OR ON THE ORDER OF, A PHYSICIAN.

VIII. INSTALLATION AND CONNECTION

A. INSTALLATION RECOMMENDATIONS

Following is a list of general recommendations when installing the EVOLUTION 3E console:

- Make sure the console is installed no more than 2 meters away from user
- The console must remain in front of user during operation to visualize every indication (value, warning)
- User must be able to hear signals from EVOLUTION 3E console properly
- User must be able to disconnect suction from aspiration tubing connector
- User must be cautious about EVOLUTION 3E console displacement due to its height, dimensions and connections
- The EVOLUTION 3E console must not be installed above or close to a radiating device.
- All applied parts must be used in normal condition of use, and on a routine surgical basis with respect to their intended use according to the instructions for use.

WARNING : The microkeratomes electric motors and turbine can reach a temperature over 41°C. Moria recommends to use the motors and rings not longer than 1 min and to prevent any contact of motors with patient tissues.

B. ELECTRICAL SETTING

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Using the screw connector, plug in the footpedal connector (#19361, 19381, 19381C, 19281J) to the rear panel of the console	Finger tight only	B8 C
2	Connect the power supply cord to the main power supply	NEVER OPERATE WITHOUT MAIN POWER SUPPLY. Before using the EVOLUTION 3E, check if position of operating voltage (115 or 230V) is appropriate to your country's power supply: - for a power supply of 100-120V: select position «115V» - for a power supply of 220-240V: select position «230V». To select position of operating voltage, see chapter XI.A. To change fuses, see chapter XI.C. If the unit has not been operated for four weeks or more, the battery level should be checked and recharged prior to use. To avoid electrical shock, the console must only be connected to a grounded power system. Position the console so that displays, connections and buttons remain visible and accessible. Always allow an easy access to the electric supply and power socket. In cases of emergency, disconnect the power socket or the supply cord connector.	B11
3	Switch on the control unit turning B9 button to "ON"	Visually check all LEDs are blinking (A3, A4, A5, A7, A9, A10, A11, A12, A14, A15, A16, A17, A18). Visually check A9 and A11 LEDs are ON. The display indicates the atmospheric pressure (approximately 760mmHg at sea level).	B9 LEDS (A9 and A11) ON



C. BATTERY TEST

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Switch EVOLUTION3E console off by turning button (B9) on «OFF» Unplug the supply cord Switch EVOLUTION3E console ON by turning button (B9) on «ON»	Check that battery charge level indicator (green LED) is switched on (wait 15 seconds)	B9 A11
2	Switch EVOLUTION 3E console by turning button (B9) on «OFF»	Check that residual gas drain switch is working (audible sound)	B9 B4
3	Plug again the supply cord to the unit	Perform surgery on main supply only.	

D. CONNECTING THE ASPIRATION TUBING

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Check the integrity of the aspiration tubing (#19138)	Visually check tubing integrity before and after set up. MORIA strongly advises against the re-use of aspiration tubing as it may damage connector, resulting in vacuum loss or damage the vacuum pumps.	
2	Connect the sterile, single use aspiration tubing (#19138) to the unit: - place the chamber in position (H1) - then get it down (H2)	Change tubing for each patient. Poorly connected tubing, either on the console or on the ring, may also cause insufficient aspiration. It is very important to check the integrity of all connections. Warning: - do not reverse the chamber position of the aspiration tubing (H3) - if liquid appears in the chamber, replace the chamber immediately (H4) - visually check tubing integrity before and after set up (H5) If moisture seeps into the console, activate vacuum to remove it. If water is sucked up, it is recommended to let the pumps run to let residual water out.	A24 H

E. VACUUM TEST

Steps	What to do	 IMPORTANT WARNINGS	Related picture																		
1	Clamp the aspiration tubing at its end (close to the ring connector)	Before connecting the suction ring on the aspiration tubing, check that internal air lumen of the ring are not obstructed by surgical residues. This obstruction, resulting from insufficient or inappropriate cleaning, would create pseudo-aspiration and risk the loss of vacuum.																			
2	Press the test button (A8) The unit will check each pump: P1, then P2. Vacuum test ends with a short audible “beep”.	This test procedure will check the efficiency of the pumps and the pressure (if gas turbine is selected). The test should be made : • prior to surgery, • whenever tubing is changed, • whenever a red indicator light is illuminated The unit can only be operated if the test is passed. The display (A7) indicates: • “PE” • then atmospheric pressure (mmHg) • then “P1” • then the pressure obtained by pump P1 (mmHg) • then “P2” • then the pressure obtained by pump P2 (mmHg) • then atmospheric pressure (mmHg) When atmospheric pressure “PE” is indicated on the display, one can hear a single “beep”. Stopping suction is indicated by a “beep” sound”.	A7 A8																		
3	The test is passed if LED (A3) is ON	The vacuum level provided by the console depends on the local atmospheric pressure, and thus on the elevation. The system displays the actual vacuum value (in mmHg). The higher the elevation, the less efficient the suction. The atmospheric pressure decreases with elevation (approx: 8 mmHg/100 m). The suction capacity is proportional to the difference in pressure between the lowest level of vacuum pump and the local atmospheric pressure. The higher the vacuum power, the smaller the value displayed. When the ring is fixated on the eye the vacuum value must be less than 250mmHg. If the value displayed is above, the unit must be immediately turned off and sent back to Moria for a servicing.  At sea level (Elevation = 0 m) Atmospheric pressure = 760 mmHg Elevation (for example 3000 m) Atmospheric pressure = 526 mmHg <table><thead><tr><th>Elevation</th><th>Atmospheric pressure mmHg</th><th>Indicative Vacuum level provided by the unit</th></tr></thead><tbody><tr><td>0 m</td><td>760</td><td>152 mmHg</td></tr><tr><td>1000 m</td><td>675</td><td>135 mmHg</td></tr><tr><td>2000 m</td><td>596</td><td>119 mmHg</td></tr><tr><td>2500 m</td><td>560</td><td>112 mmHg</td></tr><tr><td>3000 m</td><td>526</td><td>105 mmHg</td></tr></tbody></table>	Elevation	Atmospheric pressure mmHg	Indicative Vacuum level provided by the unit	0 m	760	152 mmHg	1000 m	675	135 mmHg	2000 m	596	119 mmHg	2500 m	560	112 mmHg	3000 m	526	105 mmHg	A3 ON
Elevation	Atmospheric pressure mmHg	Indicative Vacuum level provided by the unit																			
0 m	760	152 mmHg																			
1000 m	675	135 mmHg																			
2000 m	596	119 mmHg																			
2500 m	560	112 mmHg																			
3000 m	526	105 mmHg																			
4	Stop clamping the aspiration tubing. Control unit is ready to use.																				

F. LOW VACUUM FUNCTION

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Press on the vacuum footswitch (C1) then the "Low Vacuum" button (A6)	When "Low Vacuum" function is activated, LED A5 is ON. The "Low Vacuum" allows for a low level of vacuum of the suction ring on the eye. For safety reasons, when the low vacuum function is activated, the electrical motor can only be activated for a backward pass.	C1 A5 ON A6

G. MANUAL MODE (CB, CBm, CBSU, ONE, ONE USE)

Warning: Do not connect more than one keratome at a time to the EVOLUTION 3E console.

1. PNEUMATIC TURBINE CONNECTIONS

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Connect the non-copper end of the hose to the gas supply. Ensure that there is no pressure into the hose before connecting it to the console.	MORIA supplies a black hose (black is the international code for nitrogen). All questions regarding gas, gas bottles, or regulators should be addressed to the gas supplier. The hose is delivered with a copper quick connection on one end, and free on the other end. The free end has to be connected to your own gas supply system (dry and filtered medical grade air or nitrogen bottle with its dual stage gas regulator and male nozzle). MORIA recommends using: - a high-pressure dual stage regulator - dry medical grade nitrogen gas, or dry and filtered medical grade air to operate pneumatic turbines, available from your medical gas supplier MORIA does not supply the gas.	
2	Connect the Nitrogen/ Air Hose (#19120) to the control unit		B2, B3 D1, D2
3	Before opening the dry medical grade nitrogen gas tank, or dry and filtered medical grade air tank, always set the pressure to the minimum value by turning the knurled blue knob (B1) counterclockwise.	To operate the knurled blue knob (B1), pull it toward the back.	B1
4	Pressurize the system The pressure to the generator should be between 4 and 6 bars (58 and 87 psi). Incoming pressure must not exceed 14 bars.	If the pressure going IN the console exceeds 8 bars (116 psi), the console detects excess pressure and opens the high pressure valve. The dry medical grade nitrogen gas, or dry and filtered medical grade air is ejected from the turbine connector (front panel) in order to protect the internal components. If so, have the system checked by MORIA.	
5	Set the pressure in the console by turning the knob (B1) on the rear panel to obtain a pressure of 3.2 ± 0.2 bars (46 ± 3 psi)	Once the dry medical grade nitrogen gas tank, or dry and filtered medical grade air tank is opened, increase the pressure to the appropriate amount (3.2 bars): - turn the knurled blue knob (B1) clockwise to increase pressure - turn the knurled blue knob (B1) counterclockwise to decrease pressure. Once the correct pressure has been selected, lock the knob (B1) by pushing it. Do not use the unit if LED (A16) or (A18) are ON.	B1
6	Select manual mode (A13)	Press « Mode » button (A13) until green LED (A12) is ON (« Manual Mode »)	A12 A13 A15
7	Proceed a vacuum test	Refer to chapter VIII.C	A8
8	Connect the hose (#19353) to turbine	Finger tight only. Check that the hose and the connector are not damaged and that the gasket is present and undamaged on the front panel of the console and on the rear of the turbine.	A23
9	Control unit is ready to use.	Refer to user manual : - CB, CBm, CBSU for LASIK: #65005 - CBm, CBSU for artificial chamber : #65078 - ONE, One Use : #65009 Always check intraocular pressure with a tonometer (#19042) prior to using any keratome.	

2. DISASSEMBLY OF THE NITROGEN/AIR SUPPLY

Before disassembly the nitrogen/air supply, no nitrogen/air pressure must remain in the console.

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Close the gas tank.	After each surgical session, the console must be disconnected from the dry medical grade nitrogen gas supply, or dry and filtered medical grade air supply (internal circuit or tank). For safety reasons, no nitrogen or air pressure must remain in the console.	
2	Operate the turbine by pressing the Forward/ Turbine footpedal for few seconds until the turbine runs out of gas and comes to a stop.		C2
3	After surgery, switch off the vacuum by pressing the vacuum footpedal and remove the vacuum ring from the eye. Disconnect the tubing from the ring and the disposable head/ blade. Discard all in an appropriate container.		C1
4	After use, switch off the control unit and disconnect the turbine. Decontaminate and sterilize (by autoclave) the turbine.		B9
5	Press the Residual Gas Drain Switch until there is no pressure left in the console		B4
6	Disconnect the nitrogen/ air gas supply (internal circuit or tank) from the control unit.	Disconnect the gas supply hose: • first push the blue ring toward the unit (step 1), • then pull the hose towards the back (step 2)	B2–B3 D3a-D3b

H. AUTOMATED MODE (M2, M2SU, ONE USE-PLUS, EPI-K™)

Warning: Do not connect more than one keratome at a time to the EVOLUTION 3E console.

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Select the automated mode	Press « Mode » button (A13) until green LED (A14) is ON according to the selected microkeratome. « OP » appears on the display.	A13 A14 A15
2	Proceed a vacuum test	Refer to chapter VIII.C	A8
3	Plug the motor system into the control unit	To plug One Use-Plus handpieces: A19 To plug M2-M2SU handpiece: A21 To plug Epi-K™ handpiece: A22	A19 A21 A22
4	Select the advance speed (for M2, M2SU and One Use-Plus)	The advance speed 1 is slower than the advance speed 2. The blade oscillation rate is not modified. See applicable nomograms: • M2-M2SU : #65023 • One Use-Plus with plastic ring : #65041 • One Use-Plus with metal ring: #65067, 65077EN • One Use-Plus for artificial chamber : #65101	A20
5	Set the vacuum ring on the eye and press the vacuum footpedal. Assemble the keratome on the ring and activate forward footpedal to perform the cut. After the cut is performed, press backward footpedal.	Refer to user manual: • M2-M2SU: #65013 • One Use-Plus with plastic ring: #65039 • One Use-Plus with metal ring: #65040 • One Use-Plus for artificial chamber: #65106 • Epi-K™ with metal ring: #65043 • Epi-K™ with plastic ring: #65044 Always check intraocular pressure with a tonometer (#19042) prior to using any keratome.	
6	After surgery, switch off the vacuum by pressing the vacuum footpedal and remove the vacuum ring from the eye. Disconnect the tubing from the ring and the disposable head/ blade. Discard all in an appropriate container.		C1
7	After use, switch off the control unit and disconnect the electrical motor. Decontaminate the motor (see instructions in each user manual).		B9

IX. USE OF FOOTPEDAL

A. COMPATIBILITY FOOTPEDAL-CONSOLE

- Footpedals (#19361 and 19381) with serial number inferior to 7000 are not compatible with EVOLUTION 3E consoles (#19380) having a serial number superior to 7000.
- EVOLUTION 3E consoles (#19380) with serial number superior to 7000 are compatible only with footpedals (#19361, 19381, 19381C et 19381J) having a serial number superior to 7000.

Footpedal (#19361, 19381, 19381C, 19381J)	EVOLUTION 3E console (#19380)	
	5000 < sn < 7000	sn > 7000
5000 < sn	compatible	not compatible
5000 < sn < 7000	compatible	not compatible
sn > 7000	not compatible	compatible

B. COMPATIBILITY FOOTPEDAL-HANDPIECE

- Footpedal #19361 can be used with M2 and One Use-Plus handpieces, as well with turbines for ONE and CB. Footpedal #19361 can not be used with the Epi-K™ handpiece.
- Footpedal #19381 can be used with M2, One Use-Plus and Epi-K™ handpieces, as well with turbines for ONE and CB.

C. MANUAL MODE

What to do	 IMPORTANT WARNINGS	Related picture
Press the vacuum footswitch (C1) once: vacuum is ON Press again the vacuum footswitch (C1): vacuum is OFF	The system displays the actual vacuum value in mmHg (A7). When a suction ring is fixated on the eye, the vacuum value displayed must be less than 250 mmHg. If the value displayed is above 250 mmHg, please refer to chapter X.	C1
Press the forward footswitch (C2): blade oscillation	In case vacuum is not activated, the forward footswitch does not work. Blade oscillation stops once forward footswitch is not pressed any longer.	C2
Press the backward footswitch (C3): no action	Do not press footpedal (C2) and (C3) at the same time.	C3

D. AUTOMATED MODE

What to do	 IMPORTANT WARNINGS	Related picture
Press the vacuum footswitch (C1) once: vacuum is ON Press again the vacuum footswitch (C1): vacuum is OFF	The system displays the actual vacuum value in mmHg (A7). When a suction ring is fixated on the eye, the vacuum value displayed must be less than 250 mmHg. If the value displayed is above 250 mmHg, please refer to chapter X.	C1
Press the forward footswitch (C2): - blade oscillation - forward pass of microkeratome	In case vacuum is not activated, the forward footswitch does not work. When the forward footswitch is not pressed any longer : - microkeratome stops - blade oscillation stops. M2/One Use-Plus keratomes and turbines for ONE and CB: footpedals #19361 and 19381 operate at only one forward speed. Epi-K™ only (footpedal #19381): 3 forwards speeds available when press the footpedal (slow – medium – fast): 1st step: slow forward speed, with a continuous tone 2nd step: medium forward speed, with a louder tone than with 1st step 3rd step: fast forward speed, with a louder tone than with 2nd step Use B6 button to vary tone level.	C2 B6
Press the backward footswitch (C3): - backward pass of microkeratome - no blade oscillation	Pressing the backward footswitch induces a backward pass of the microkeratome even if: - vacuum is OFF "Low Vacuum" option is ON Do not press footpedal (C2) and (C3) at the same time.	C3

X. TROUBLE SHOOTING

A. POWER SUPPLY

Message	Reason	Intervention
LED (A9-A11) ON + Audible tone	Loss of main power supply during surgery AND Battery charged (>90%) Control unit runs on battery supply	The system runs automatically on the backup battery: check LED A11 is OFF. At the end of the surgery case, connect the console to the main power supply to recharge the battery immediately: check LED A11 is ON.
LED (A10) ON + Audible tone	Loss of Main Power supply: Main power supply is defective during surgery AND Battery discharged Control unit runs on a discharged battery supply	Connect the console to the main power supply to recharge the battery immediately.

B. VACUUM LOSS

Message	Reason	Intervention
LED (A3–A4) ON + audible tone	Loss of vacuum has been detected on pump P1 or P2. Display (A7) shows « P1 » or « P2 » depending on the defective pump, then current atmospheric value.	One of the two pumps has failed. The second pump will automatically and instantaneously take over to provide vacuum power. Even if there is a vacuum loss, the microkeratome will continue to function. At the end of the case, run a test. Then if problem still persists, check that aspiration tubing is well clamped, then run a test. Then if problem still persists, check aspiration tubing and change it if any doubt, then run a test. If unable to solve the problem, do not use the console and contact your distributor or MORIA to have it checked.
LED (A4) ON + double audible tone	Loss of vacuum has been detected during surgery on pump P1 and P2. Displays (A7) shows simultaneously P1/P2 then "___"	The two pumps have failed. Then if problem still persists, check that aspiration tubing is well clamped, then run a test. Then if problem still persists, check aspiration tubing and change it if any doubt, then run a test. If unable to solve the problem, do not use the console and contact your distributor or MORIA to have it checked. Note: Even if there is a vacuum loss during the case, the microkeratome will continue to function.

C. RESIDUAL VACUUM

If the vacuum cannot be interrupted by pressing the vacuum footpedal while suction ring is on patient's eye, cut or disconnect the vacuum tubing from the console to decrease pressure.
Then lift keratome carefully taking care not to harm the corneal hinge.

D. TURBINE PRESSURE (MANUAL MODE)

Message	Reason	Intervention
LED (A16) ON + audible tone	Loss of turbine pressure	The yellow LED will turn ON (as soon as the pressure is below 2.5 bars). At the end of the case, check the system and run a test procedure. Note: The loss of turbine pressure will decrease the blade oscillation rate.
Nitrogen/Air going out by the hose	Too much turbine pressure	If the pressure going in the console exceeds 8 bars (116 psi), it may damage the unit. MORIA strongly recommends sending back the console for a complete check-up. The pressure sensor board and other components may have been damaged and may alter the performance of the turbines and/or alter the information on the display panel.

XI. CARE & MAINTENANCE

In cases of unusual vibrations or noises, instability, unusual performance, or failure to pass the vacuum test, do not use the unit and contact MORIA or your distributor for servicing.

MORIA strongly recommends having the system thoroughly inspected by MORIA on a routine basis every year.

WARNING:

Wipe the control unit with a lint-free surgical wipe moistened with disinfecting solutions..

Do not use solvents, abrasives, fluids, or solutions on the control unit and footpedals.

Do not pour liquids or solutions on the control unit, footpedals and connectors.

Please contact MORIA for any other servicing.

A. HOW TO SELECT THE OPERATING VOLTAGE (115V OR 230V POSITION)

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Switch OFF the control unit.	Before opening the control unit, switch OFF the control unit and remove the cord.	B9
2	Gently open the door.		B10 E1-E2-E3
3	Remove the selector (A) and select the right voltage	For a power supply of 100-120V: select position «115V». For a power supply of 220-240V: select position «230V».	E4
4	Reposition the selector and close the door		E5

B. INSTALLING A NEW BATTERY

If the battery can not be charged, or does not hold a charge, a new battery must be installed in the console.
See draws F.

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Switch OFF the control unit.	Before opening the control unit, switch off the control unit and remove the cord.	B9
2	Completely unscrew the 2 screws located at the back of the control unit		F1 – B14
3	Slide the cover toward the front		F2
4	Lift the cover from back to front in 2 steps		F3
5	Unplug the battery		F4
6	Unscrew the 2 screws		F5
7	Remove the metal piece (while keeping screws in place)		F6
8	Lift the battery with the plastic strap	Maintenance guide requires returning out of order batteries to Moria for treatment.	F7
9	Replace with a new battery, reposition and screw in the metal piece, plug in the new battery, replace and screw down the cover	Battery: 12 V - 7 Ah (Pb) Use only the original manufacturer's parts. The use of a different battery may result in severe damage to the unit, or in malfunction. To recycle the battery, please contact your dealer or MORIA. NEVER USE THE CONTROL UNIT WITHOUT A BATTERY INSTALLED INSIDE AND PLUGGED IN. Important: the used battery is classified as waste, which must be dealt with according to the specific regulations for collection, treatment, recycling and disposal set out in directive 2006/66/EC of 9 September 2006 relating to batteries and accumulators and waste batteries and accumulators. As such the battery must be disposed of at an appropriate collection point.	F8
10	Test battery	Refer to chapter VIII.C	

C. REPLACING FUSES

Steps	What to do	 IMPORTANT WARNINGS	Related picture
1	Switch OFF the control unit.	Before opening the control unit, switch OFF the control unit and remove the cord.	B9
2	Gently open the door with screwdriver		B10 G1-G2
3	Remove the cartridge and change the 2 fuses	500 mA high switching power	G3-G4-G5
4	Close the door		G1

XII. WARRANTY

A. SCOPE OF WARRANTY

DESCRIPTION	MORIA reference
EVOLUTION 3E Console (Except Battery)	19380
EVOLUTION 3-3E Control Footpedal	19361
EVOLUTION 3E Control Footpedal Epi-K™	19381
EVOLUTION 3E Footswitch - China	19381C
EVOLUTION 3E footswitch - Japan	19381J

- The above items as well as spare parts and labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, after having previously been decontaminated.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.

B. NON-APPLICATION OF WARRANTY

The warranty will not be applicable under any of the following conditions:

- Defects or malfunction that occur out of the warranty period (paragraph XII.C.).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- The use of supplies, spare parts, or accessories other than those supplied by MORIA (for example: blades and tubing not supplied by MORIA).
- Disassembly, modification, or intervention carried out on the devices by a person not authorized by MORIA.

C. WARRANTY PERIOD

- The warranty takes effect on the date the material is dispatched.
- The duration of the warranty is 12 months from the date of effect.

D. LIABILITY

- The liability of MORIA is limited to the supply of the services mentioned in paragraph XII.A. MORIA will not be held responsible for any direct or indirect damage notably financial, suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre, France will have sole jurisdiction.

XIII. DRAWINGS

A. CONTROL UNIT (FRONT)

Figure	Description
1	Cover
2	Chamber holder
3	Pump indicator – green LED
4	Pump indicator – red LED
5	Low Vacuum LED indicator
6	Low Vacuum switch
7	Vacuum Level Indicator
8	“TEST” switch
9	Battery Level Indicators : Acceptable – green LED
10	Battery Level Indicators: Low – red LED
11	Main Power Supply Indicator – green LED
12	Manual mode indicator
13	Mode selection switch
14	Mode Display : automated
15	Mode Display and pressure display
16	Low pressure display – yellow LED
17	Correct pressure display – green LED
18	High pressure display – yellow LED
19	One Use-Plus blue connector
20	Speed 1/Speed 2 Switch (for M2, M2SU and One Use-Plus)
21	M2 grey connector
22	Epi-K™ yellow connector
23	Turbine connector
24	Aspiration tubing connector

B. CONTROL UNIT (BACK)

Figure	Description
1	Control knob to adjust gas pressure to the turbine
2	Quick Connector for Nitrogen/Air supply
3	Connector release
4	Residual Gas Drain Switch
5	MORIA reference and Serial Number
6	Volume of Epi-K™ Footpedal Steps
7	CE mark, manufacturer address
8	Footpedal Connector
9	Main OFF (0) /ON (1) Switch
10	Fuse Housing and 115V – 230V Selector
11	Main Power Supply Connector
12	Power supply informations
13	Pumps air release
14	Screws

C. FOOTPEDAL

Figure	Description
1	Vacuum footpedal
2	Forward keratome pass or Turbine ON
3	Backward pass
4	IPX6 (serial number above 8000)

D. CONNECTING AND DISCONNECTING NITROGEN/ AIR HOSE (#19120)

Figure	Description
1	To connect : approach the hose...
2	... then screw it
3a	To disconnect : first push the blue ring toward the unit...
3b	... then pull the hose towards the back

E. HOW TO SELECT THE OPERATING VOLTAGE (115V OR 230V POSITION)

Figure	Description
1	Switch OFF the control unit.
2	Gently open the door.
3-4	Remove the selector and select the right voltage
5	Reposition the selector and close the door

F. CHANGING THE BATTERY

Figure	Description
1	Unscrew completely the 2 screws located at the back of the control unit
2	Slide the cover towards the front
3	Lift the cover from back to front in 2 steps (1 – 2)
4	Unplug the battery
5	Unscrew the 2 screws
6	Remove the metal piece (while keeping screws in place)
7	Lift the battery with the plastic strap
8	Replace with a new battery, reposition and screw the metal piece, plug the new battery, replace and screw the cover
9	Discard the battery in an appropriate collection bin. The used battery is classified as waste, which must be dealt with according to the specific regulations for collection, treatment, recycling and disposal set out in directive 2006/66/EC of 9 September 2006 relating to batteries and accumulators and waste batteries and accumulators. As such the battery must be disposed of at an appropriate collection point.

G. CHANGING FUSES

Figure	Description
1	Switch OFF the control unit.
2	Gently open the door with screwdriver
3	Remove the cartridge
4	Change the 2 fuses
5	Replace the 2 fuses

H. INSTALLING TUBING

Figure	Description
1	Place the chamber in position
2	Get the chamber down in the chamber maintainer
3	Do not reverse chamber position of the aspiration tubing
4	If liquid appears in the chamber, replace it immediately
5	Visually check tubing integrity before and after set up

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MORIA S.A. – 15, rue Georges Besse - 92160 ANTONY – France