#65065

EN

Instructions for use of the One® Single Use Corneal Vacuum Punch #17200Dxxx

Indication for use

The Moria ONE® Single-Use Corneal Vacuum Punch is designed for performing corneal grafts (lamellar or penetrating keratoplasty).

Device components

- Cutting block base (A) (see Fig. 1)
- Trephine blade housing (B) (see Fig. 1)
- Syringe (C) (see Fig. 1)

The most recent version of this document is available on the Moria website (<http://www.moria-surgical.com>) via the restricted access heading «Documentation/Library». Translations are also available in other languages under that same heading. Registration required.


Warning

- These instructions explain only how to use this medical device and do not describe the steps necessary to perform a corneal graft. The surgical procedure remains the responsibility of the surgeon experienced in performing corneal grafts.
- Check the integrity of the pouch.
- Check also the expiration date on the packaging.
- The ONE® single-use corneal vacuum punch is sterile and ready for use. Do not use it if the pouch in which it is packaged is damaged.
- Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performances and exposes the patient to adverse events.

Directions

- 1 / Check the blade diameter which is stamped on top of the trephine blade housing.
- 2 / Remove the sterile blister from the pouch.
- 3 / Remove the cutting block base, trephine blade housing and the syringe from the blister.
- 4 / Connect the tubing to the syringe with the Luer™ lock connection (see Fig 2).
- 5 / Check the integrity of the trephine blade under the microscope.
- 6 / Press the piston of the syringe half way and hold.
- 7 / Place the graft on the cutting block base, endothelial side up.
- 8 / Center the graft using the 4 cardinal holes diametrically opposite of 9.5 mm.
- 9 / Release the piston of the syringe to create the vacuum that will hold the graft in position on the cutting block base.
Check that the graft is still well centered.
If the graft is not well centered at the end of this preparation, repeat the steps 6, 7, 8 and 9 to obtain a perfect centering.
- 10 / Place the trephine blade housing on the cutting block base, let it lower down by gravity to avoid setting it with a wrong angle (risk of blockage) (see Fig 3) and press down to cut the graft (see Fig 4). While holding the assembled punch, grasp the sclera, through the "top lateral" window, with a tooth forceps and turn it gently to finalize the cut (see Fig 5).
- 11 / Pull gently upward the trephine housing blade from the cutting block base. The donor button should be kept in place on the cutting block base thanks to the aspiration. Then remove completely the trephine housing blade.
- 12 / Disassemble the syringe from the tubing to release the vacuum. The donor button is now ready to be used on the patient.
- 13 / At the end of the surgical procedure, discard the ONE® single-use vacuum punch into an appropriate container for destruction (according to the regulations in force).

Explanation of symbols (General list of symbols used by Moria)

 XXXXXX	CATALOGUE REFERENCE
	QUANTITY
 YYYY/MM/DD	USE BY
 XXXXXX	BATCH CODE
	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
 XXXXXXX	MANUFACTURER
 YYYY/MM/DD	DATE OF MANUFACTURE
 XXXXXX	CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)
 XXXXXX	CONSULT OPERATING INSTRUCTIONS
	DO NOT USE IF PACKAGE IS DAMAGED
	STORAGE CONDITIONS: • TEMPERATURE: 15°C – 35°C / 59°F – 95°F • HUMIDITY RATE: 30% – 70%
	KEEP DRY
	DISCARD IN A SEPARATE COLLECTION BIN
	CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.

List of references

Reference: # 17200Dxxx	
17200D600	Single Use Corneal Vacuum Punch D600
17200D650	Single Use Corneal Vacuum Punch D650
17200D675	Single Use Corneal Vacuum Punch D675
17200D700	Single Use Corneal Vacuum Punch D700
17200D725	Single Use Corneal Vacuum Punch D725
17200D750	Single Use Corneal Vacuum Punch D750
17200D775	Single Use Corneal Vacuum Punch D775
17200D800	Single Use Corneal Vacuum Punch D800
17200D825	Single Use Corneal Vacuum Punch D825
17200D850	Single Use Corneal Vacuum Punch D850
17200D875	Single Use Corneal Vacuum Punch D875
17200D900	Single Use Corneal Vacuum Punch D900
17200D950	Single Use Corneal Vacuum Punch D950
Packaging	1 per box, sterile, single-use

Accessories

Reference	17202Dxxx
Description	One® - Adjustable Vaccum Trephine Diameters: 6.50 to 9.00 mm with 0.25 mm increments. Also 6.00, 9.50 and 10.00 mm diameters
Packaging	1 per box, sterile, single-use

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Referência: #17200Dxxx	
17200D600	ONE® Single Use Corneal Vacuum Punch D600
17200D650	ONE® Single Use Corneal Vacuum Punch D650
17200D675	ONE® Single Use Corneal Vacuum Punch D675
17200D700	ONE® Single Use Corneal Vacuum Punch D700
17200D725	ONE® Single Use Corneal Vacuum Punch D725
17200D750	ONE® Single Use Corneal Vacuum Punch D750
17200D775	ONE® Single Use Corneal Vacuum Punch D775
17200D800	ONE® Single Use Corneal Vacuum Punch D800
17200D825	ONE® Single Use Corneal Vacuum Punch D825
17200D850	ONE® Single Use Corneal Vacuum Punch D850
17200D875	ONE® Single Use Corneal Vacuum Punch D875
17200D900	ONE® Single Use Corneal Vacuum Punch D900
17200D950	ONE® Single Use Corneal Vacuum Punch D950
Embalagem	1 por caixa, esterilizado, utilização única

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MORIA Inc